



DASH Lite

YOUR ULTIMATE TRAVEL COMPANION

BEST VALUE FOR MONEY

THE PERFECT MATCH OF
COMFORT & WEIGHT



PLUSH
SEATING



WEIGHS ONLY
46 LBS


COMFORTABLE RIDING EXPERIENCE

No need to sacrifice comfort for lightweight, have them both with DASH. Featuring ample hip room, a supportive backrest, thick plush seating, soft armrests, and large shock absorbing rear wheels, DASH guarantees a comfortable riding experience.



FOLD IT, LOAD IT, AND GO!

Boasting a simple one-step fold sequence, DASH's unique folded design allows it to be rolled like a suitcase & fits in the back of any vehicle for convenient travel. DASH is lightweight and compact allowing many users to lift it on their own without the need for special transport vehicles, vans, or expensive lifts. DASH is the ultimate travel power chair!



Intelligently Designed

KEY SPECIFICATIONS



WEIGHT CAPACITY
265 LBS



DRIVING RANGE
9.5 MI / 15 KM



BATTERY TYPE
AIRLINE APPROVED



INCLINE CAPABILITY
9° / 15.9%



FOLDED DIMENSIONS
30" L X 24" W X 13.5" H



UNFOLDED DIMENSIONS
38" L X 24" W X 37" H



SEATING WIDTH
19"



BACK REST HEIGHT
19"

When you choose Travel Buggy, you're not just getting a mobility solution; you're becoming part of a vibrant community of adventurers, explorers, and advocates for accessible travel. It is time to **SEE YOUR WORLD AGAIN** with Travel Buggy.



TRAVEL BUGGY WAS FOUNDED TO PROVIDE
TRAVEL-FRIENDLY POWER MOBILITY EQUIPMENT

We offer more accessories than anyone else. From essentials like cup holders & extra storage to vibrant seating colors. Navigate Your World in Style with DASH!



FOUR VIBRANT COLORS