

WORLD'S LIGHTEST HEAVY DUTY WHEELCHAIR!

NEW

Feather™ CHAIR HD

Where heavy duty meets light weight.



The Feather Chair HD boasts a 22" or 24" wide seat, 350 lbs. or 450 lbs. capacity and weighs only 22 lbs., **Pop off the wheels**, and it weighs only 15 lbs.! Swing-away footrests and anti-tippers offer comfort and safety. The **back folds down** for extra added convenience during storage or travel.

FEATURES:

- only 15 lbs without wheels
- 22" or 24" wide seat
- folding frame
- anti-tippers
- optional elevating leg rests


SPECIFICATIONS

Weight Capacity	350 /450 lbs.		
Product Weight	22 / 23 lbs.	Front Wheel Size	8"
Weight Without Wheels	15 /16 lbs.	Rear Wheel Size	24"
Seat Width	22" / 24"	Axle Type	Fixed, Quick Release
Seat Depth	16"	Armrest Options	Full Length
Seat-to-Floor Height	19"	Armrest Height	9"
Overall Width	32"	Back Height	16"
Overall Length	34"	Upholstery Type	Nylon
Overall Height	28"	Upholstery Color	Black w/ Green Trim or Black w/ Silver Trim
Folded Depth	13"	Footrests	Swingaway
Folded Height	28"	Warranty on Frame	Lifetime Limited
Folded Width	15"	Warranty on Wearables	Six Months Limited

Lightweight



Foldable



Locking Brakes



Arrives Assembled

