

GLIDE CARE G5 *by Comfygo Mobility*


THE LOW COST - 45 POUND - REMOTE CONTROLLED - FOLDABLE POWER CHAIR



Battery



Motor




Airline



Cruise



Net Weight



Range



Max Speed



Max Slope



Width



Capacity



10AH

2x250W

APPROVED

APPROVED

45lbs

12mi

4mph

9deg

22in

265lb



SPECIFICATION	G5 DETAIL
Overall Dimensions (L x W x H)	36"x22"x37" (920x565x940 mm)
Folded Dimensions (L x W x H)	29.5"x22"x12" (750x565x300 mm)
Seat Height (from ground)	19" (480 mm)
Armrest Height (from seat)	9" (230 mm)
Backrest Height (from seat)	19" (490 mm)
Max Loading Weight	265 lbs (120 kg)
Net Weight without Battery	46 lbs (20.8 kg)
Net Weight with Battery	50 lbs (22.6 kg)
Battery	24V 10AH lithium-ion
Battery Charge	AC 100-240V, 50Hz, 1.5-5A
Battery Charge Time	4 to 8 hours
Range with 12AH Battery	Up to 10 mile (16 km)
Airline and Cruise Approval	Approved
Max Speed	4 mph (6 km/h)
Motor	2 x 250W Brushed (Total 500W)
Controller	24V 35A
Obstacle Climbing	1.5" (400 mm)
Climbing Angle Max Degrees	10°
Max Safe Slope	9°
Static Stability	9°
Brake	Intelligent electromagnetic brake
Turning Radius	≤ 47" (1200 mm)
Reversing Width	≤ 59" (1500 mm)
Front Wheel Diameter	7" (178 mm)