

G5 Pro_{by Comfygo Mobility}

45 POUND - REMOTE CONTROLLED - FOLDABLE POWER CHAIR - SUITCASE HANDLE



Battery



10AH

Motor



2x250W

Airline



APPROVED

Cruise



APPROVED

Net Weight



45lbs

Range



13mi

Max Speed



4mph

Max Slope



9deg

Width



22in

Capacity



265lb



G5 PRO TECHNICAL SPECS

Parameters	G5 Pro
Overall Dimensions (L x W x H)	36"x22"x37" (920x565x940 mm)
Folded Dimensions (L x W x H)	29.5"x22"x12" (750x565x300 mm)
Seat Height (from ground)	19" (480 mm)
Armrest Height (from seat)	9" (230 mm)
Backrest Height (from seat)	19" (490 mm)
Max Loading Weight	265 lbs (120 kg)
Net Weight without Battery	40 lbs (18 kg)
Net Weight with Battery	45 lbs (20.5 kg)
Battery	24V 10AH lithium-ion
Battery Charge	AC 100-240V, 50Hz, 1.5-5A
Battery Charge Time	4 to 8 hours
Range with 12AH Battery	Up to 13 mile (20 km)*
Airline and Cruise Approval*	Approved**
Max Speed	4 mph (6 km/h)
Motor	2 * 250W Brushless (Total 500W)
Controller	24V 40A Brushless
Obstacle Climbing	1.5" (400 mm)
Climbing Angle Max Degrees	10°
Max Safe Slope	9°
Static Stability	9°
Brake	Intelligent electromagnetic brake
Turning Radius	≤ 47" (1200 mm)
Reversing Width	≤ 59" (1500 mm)
Front Wheel Diameter	7" (178 mm)
Rear Wheel Diameter	11" (280 mm)
Seating Area Width	18" (460 mm) ***
Cushion Size (width x depth)	17.5" (450 mm) x 16.9" (430 mm)
Controller Type	Detachable All-Direction Joystick