

GATEWAY® 3G Solid Surface Portable Ramp Technical Specifications

OVERALL

Material : 6000 Series Aluminum


Slip-resistant Surface : Extruded or SAFEGUARD Extra Fine

Direct Gritting

Weight Capacity : 1,000 lbs.

Warranty Term : Lifetime Limited*

Use : Wheelchairs & Scooters



GATEWAY 3G 3-FT RAMP

Model No. : GATEWAY3G 03

Usable Size (L x W) : 3' x 36"

Overall Size (L x W x H) : 3' x 37" x 3.5"

Curb Height : 2"
Weight** : 24 lbs.
Usage Guidelines (Rise)*** : Up to 4"

GATEWAY 3G 4-FT RAMP

Model No. : GATEWAY3G 04

Usable Size $(L \times W)$: $4' \times 36''$

Overall Size (L x W x H) : 4' x 37" x 3.5"

Curb Height : 2"
Weight** : 30 lbs.
Usage Guidelines (Rise)*** : Up to 6"

GATEWAY 3G 5-FT RAMP

Model No. : GATEWAY3G 05

Usable Size (L x W) : $5' \times 36''$

Curb Height : 2"
Weight** : 36 lbs.
Usage Guidelines (Rise)*** : Up to 7"

GATEWAY 3G 6-FT RAMP

Model No. : GATEWAY3G 06

Usable Size (L x W) : $6' \times 36''$

Overall Size (L x W x H) : 6' x 37" x 3.5"

Curb Height : 2"
Weight** : 42 lbs.
Usage Guidelines (Rise)*** : Up to 9"

* Warranty excludes powder-coated finish.

** Weights and dimensions are approximate and subject to change.

GATEWAY 3G 7-FT RAMP

Model No. : GATEWAY3G 07

Usable Size (L x W) : 7' x 36"

Overall Size (L x W x H) : 7' x 37" x 3.5"

Curb Height : 2"
Weight** : 50 lbs.
Usage Guidelines (Rise)*** : Up to 10"

GATEWAY 3G 8-FT RAMP

Model No. : GATEWAY3G 08

Usable Size (L x W) : 8' x 36"

Overall Size (L x W x H) : 8' x 37" x 3.5"

Curb Height : 2"
Weight** : 54 lbs.
Usage Guidelines (Rise)*** : Up to 12"

GATEWAY 3G 9-FT RAMP

Model No. : GATEWAY3G 09

Usable Size (L x W) : 9' x 36" Overall Size (L x W x H) : 9' x 37" x 5.5"


Curb Height : 4"
Weight** : 68 lbs.
Usage Guidelines (Rise)*** : Up to 13"

GATEWAY 3G 10-FT RAMP

Model No. : GATEWAY3G 10

Usable Size (L x W) : 10' x 36" Overall Size (L x W x H) : 10' x 37" x 5.5"

Curb Height : 4"
Weight** : 72 lbs.
Usage Guidelines (Rise)*** : Up to 15"







^{***} Consult your mobility equipment manufacturer for proper degree of incline; never exceed their recommendations.



GATEWAY® 3G Solid Surface Portable Ramp w/Two-line Handrails Technical Specifications

OVERALL

Loops : Included

Material : 6000 Series Aluminum

Slip-resistant Surface : Extruded or SAFEGUARD Extra Fine


Direct Gritting

Handrail Surface : Architectural Finish or Powder-Coated

Weight Capacity : 1,000 lbs.

Warranty Term : Lifetime Limited*

Use : Wheelchairs & Scooters



GATEWAY 3G 3-FT RAMP W/TWO-LINE HANDRAILS

Model No. : GATEWAY3G 03TL Usable Size (L x W x H) : 3' x 36" x 34.5" Overall Size (L x W x H) : 3' x 41.75" x 36"

Curb Height : 2"
Handrail Length : 2'
Weight** : 44 lbs.
Usage Guidelines (Rise)*** : Up to 4"

GATEWAY 3G 4-FT RAMP W/TWO-LINE HANDRAILS

Model No. : GATEWAY3G 04TL
Usable Size (L x W x H) : 4' x 36" x 34.5"

Overall Size (L x W x H) : 4' x 41.75" x 36"

Curb Height : 2"
Handrail Length : 3'
Weight** : 52 lbs.
Usage Guidelines (Rise)*** : Up to 6"

GATEWAY 3G 5-FT RAMP W/TWO-LINE HANDRAILS

Model No. : GATEWAY3G 05TL
Usable Size (L x W x H) : 5' x 36" x 34.5"
Overall Size (L x W x H) : 5' x 41.75" x 36"

Curb Height : 2"
Handrail Length : 4'
Weight** : 60 lbs.
Usage Guidelines (Rise)*** : Up to 7"

GATEWAY 3G 6-FT RAMP W/TWO-LINE HANDRAILS

 Model No.
 : GATEWAY3G 06TL

 Usable Size (L x W x H)
 : 6' x 36" x 34.5"

 Overall Size (L x W x H)
 : 6' x 41.75" x 36"

Curb Height : 2"
Handrail Length : 5'
Weight** : 68 lbs.
Usage Guidelines (Rise)*** : Up to 9"

GATEWAY 3G 7-FT RAMP W/TWO-LINE HANDRAILS

 Model No.
 : GATEWAY3G 07TL

 Usable Size (L x W x H)
 : 7' x 36" x 34.5"

 Overall Size (L x W x H)
 : 7' x 41.75" x 36"

Curb Height : 2"
Handrail Length : 6'
Weight** : 74 lbs.
Usage Guidelines (Rise)*** : Up to 10"

GATEWAY 3G 8-FT RAMP W/TWO-LINE HANDRAILS

 Model No.
 : GATEWAY3G 08TL

 Usable Size (L x W x H)
 : 8' x 36" x 34.5"

 Overall Size (L x W x H)
 : 8' x 41.75" x 36"

Curb Height : 2"
Handrail Length : 7'
Weight** : 82 lbs.
Usage Guidelines (Rise)*** : Up to 12"

GATEWAY 3G 9-FT RAMP W/TWO-LINE HANDRAILS

 Model No.
 : GATEWAY3G 09TL

 Usable Size (L x W x H)
 : 9' x 36" x 34.5"

 Overall Size (L x W x H)
 : 9' x 41.75" x 36"

Curb Height : 4"
Handrail Length : 8'
Weight** : 96 lbs.
Usage Guidelines (Rise)*** : Up to 13"

GATEWAY 3G 10-FT RAMP W/TWO-LINE HANDRAILS

Model No. : GATEWAY3G 10TL
Usable Size (L x W x H) : 10' x 36" x 34.5"

Overall Size (L x W x H) : 10' x 41.75" x 36"

Curb Height : 4"
Handrail Length : 9'
Weight** : 106 lbs.
Usage Guidelines (Rise)*** : Up to 15"

^{*} Warranty excludes powder-coated finish.

^{**} Weights and dimensions are approximate and subject to change.

^{***} Consult your mobility equipment manufacturer for proper degree of incline; never exceed their recommendations.



GATEWAY® 3G Solid Surface Portable Ramp w/Vertical Pickets Technical Specifications

OVERALL

Loops : Included

Material : 6000 Series Aluminum

Slip-resistant Surface : Extruded or SAFEGUARD Extra


Fine Direct Gritting

Handrail Surface : Architectural Finish or

Powder-Coated

Weight Capacity : 1,000 lbs.

Warranty Term : Lifetime Limited*
Use : Wheelchairs & Scooters



GATEWAY 3G 3-FT RAMP W/VERTICAL PICKETS

Model No. : GATEWAY3G 03VP Usable Size (L x W x H) : 3' x 36" x 34.5" Overall Size (L x W x H) : 3' x 41.75" x 36"

Curb Height : 2"
Handrail Length : 2'
Weight** : 46 lbs.
Usage Guidelines (Rise)*** : Up to 4"

GATEWAY 3G 4-FT RAMP W/VERTICAL PICKETS

Model No. : GATEWAY3G 04VP Usable Size (L x W x H) : 4' x 36" x 34.5" Overall Size (L x W x H) : 4' x 41.75" x 36"

Curb Height : 2"
Handrail Length : 3'
Weight** : 58 lbs.
Usage Guidelines (Rise)*** : Up to 6"

GATEWAY 3G 5-FT RAMP W/VERTICAL PICKETS

Model No. : GATEWAY3G 05VP Usable Size (L x W x H) : 5' x 36" x 34.5" Overall Size (L x W x H) : 5' x 41.75" x 36"

Curb Height : 2"
Handrail Length : 4'
Weight** : 68 lbs.
Usage Guidelines (Rise)*** : Up to 7"

GATEWAY 3G 6-FT RAMP W/VERTICAL PICKETS

Model No. : GATEWAY3G 06VP Usable Size (L x W x H) : 6' x 36" x 34.5" Overall Size (L x W x H) : 6' x 41.75" x 36"

Curb Height : 2"
Handrail Length : 5'
Weight** : 78 lbs.
Usage Guidelines (Rise)*** : Up to 9"

GATEWAY 3G 7-FT RAMP W/VERTICAL PICKETS

Model No. : GATEWAY3G 07VP Usable Size (L x W x H) : 7' x 36" x 34.5" Overall Size (L x W x H) : 7' x 41.75" x 36"

Curb Height : 2"
Handrail Length : 6'
Weight** : 90 lbs.
Usage Guidelines (Rise)*** : Up to 10"

GATEWAY 3G 8-FT RAMP W/VERTICAL PICKETS

Curb Height : 2"
Handrail Length : 7'
Weight** : 98 lbs.
Usage Guidelines (Rise)*** : Up to 12"

GATEWAY 3G 9-FT RAMP W/VERTICAL PICKETS

Model No. : GATEWAY3G 09VP Usable Size (L x W x H) : 9' x 36" x 34.5" Overall Size (L x W x H) : 9' x 41.75" x 36"

Curb Height : 4"
Handrail Length : 8'
Weight** : 118 lbs.
Usage Guidelines (Rise)*** : Up to 13"

GATEWAY 3G 10-FT RAMP W/VERTICAL PICKETS

Model No. : GATEWAY3G 10VP Usable Size (L x W x H) : 10' x 36" x 34.5" Overall Size (L x W x H) : 10' x 41.75" x 36"

Curb Height : 4"
Handrail Length : 9'
Weight** : 128 lbs.
Usage Guidelines (Rise)*** : Up to 15"

^{*} Warranty excludes powder-coated finish.

^{**} Weights and dimensions are approximate and subject to change.

^{***} Consult your mobility equipment manufacturer for proper degree of incline; never exceed their recommendations.