

# AL500 / AL500HD

## UNIVERSAL POWER WHEELCHAIR LIFT



### DESIGNED FOR POWER WHEELCHAIRS

The AL500 series lifts are among the most popular outside platform lifts. Featuring Harmar's Stars -N- Stripes aluminum deck design, dedicated to our veterans and commemorating Harmar's heritage as an American company.

The AL500 and AL500HD are easy to operate and built to last for years of trouble-free use. The universal design safely transports the power chair without making chair modifications.

- Easy, one-switch operation
- License Plate Light
- Key Switch and Courtesy Light
- 4- self-tensioning Q'Straint Retractors, no chair modifications required
- Manual Crank backup and License Plate Mount included
- Swing-Away option available

### MULTIVIEW BACKUP CAMERA

The backup camera that won't be blocked by lift platforms or mobility devices!



Camera mounts to the edge of the platform



In-vehicle moveable display monitor or rearview mirror



Camera lenses are fixed in two positions to automatically provide a clear view when platform is stowed or deployed with a mobility device

*Warranty is 1-year (non-transferable)*



*I am a caretaker for my father and this lift made it possible for me to take my dad to go out to eat and meet with friends and family. He now feels a sense of freedom and enjoys trips big and small.*

## AL500 PRODUCT SPECIFICATIONS

<b>LIFTING CAPACITY</b>	350 lb
<b>ACCOMMODATES</b>	Power wheelchairs
<b>INSTALLED WEIGHT</b>	75 lb
<b>PLATFORM SIZE</b>	26" x 39"
<b>HITCH CLASS</b>	II or greater
<b>POWERED</b>	Fully powered
<b>WARRANTY</b>	3-year parts
<b>INCLUDED</b>	<ul style="list-style-type: none"> <li>- AL107 hitch post</li> <li>- AL123 class III hitch adapter</li> <li>- License plate mount</li> </ul>
<b>SECUREMENT</b>	4 self-tensioning Q'Straint retractors
<b>AVAILABLE OPTIONS</b>	<ul style="list-style-type: none"> <li>- Battery pack</li> <li>- Cover</li> <li>- License plate riser</li> <li>- MultiView Backup Camera</li> <li>- Swing-away</li> </ul>

## AL500HD PRODUCT SPECIFICATIONS

<b>LIFTING CAPACITY</b>	400 lb
<b>ACCOMMODATES</b>	Power wheelchairs
<b>INSTALLED WEIGHT</b>	77 lb
<b>PLATFORM SIZE</b>	26" x 39"
<b>HITCH CLASS</b>	III or greater
<b>POWERED</b>	Fully powered
<b>WARRANTY</b>	3-year parts
<b>INCLUDED</b>	<ul style="list-style-type: none"> <li>- AL107 hitch post</li> <li>- AL123HD class III heavy duty hitch adapter</li> <li>- License plate mount</li> </ul>
<b>SECUREMENT</b>	4 self-tensioning Q'Straint retractors
<b>AVAILABLE OPTIONS</b>	<ul style="list-style-type: none"> <li>- Battery pack</li> <li>- Cover</li> <li>- License plate riser</li> <li>- MultiView Backup Camera</li> <li>- Swing-away</li> </ul>