


MAJESTIC IQ-CARBON

Remote-Controlled Carbon Fiber Lightweight Electric Wheelchair



Battery




24V - 12Ah
24V - 10Ah

Range




Up to 13.6 mi (22km)*
Up to 11 mi (18km)*

Max Speed



4 mph
(6 km/h)

Motor



2 x 250W
Brushless

Airline | Cruise



Approved
with 10Ah / 12Ah

Capacity



300 lbs
(136 kg)

LCD Display



Remote Control



*The actual range of our product can vary due to factors such as weather conditions, load weight, and terrain type. Extreme temperatures, heavy loads, and challenging terrains can all impact battery efficiency and reduce travel distance.



Premium carbon fiber frame for ultra-lightweight strength, made for travel..
Features dual 250W whisper-quiet motors, up to 13.6-mile range, airline approved,
a sleek LCD display joystick, USB charging, and headlight – built for comfort,
control, and high performance.



Black



Taba



Red



Coffee



Gray



Blue



Standard



Purple



Teal



35"
(90 cm)

20"
(51 cm)

39"
(99.5 cm)

23.5"
(60 cm)

Seating Width
21" (53 cm)

Max Climb	Net Weight	Heaviest Piece	Cushion Size	Chair Width	Turn Radius
					
10° **	45 lbs. (20.5 kg)	38.4 lbs. (17.4 kg)	Width: 17.7" (45 cm) Depth: 17.7" (45 cm)	23.5" (60 cm)	35.5" (90 cm)



**The steepest slope the chair can climb under test conditions.

makes life easier



10