

# MAJESTIC IQ-8000



**COMFYGO**  
mobility



12AH or 20AH  
lithium-ion



2 \* 250 watts



Approved \*



Approved \*



58.5 lbs



Up to 19 miles\*



4 mph




6 degrees



23.5" or 26"



330 lbs



**Seating Width\*\***  
17.5" or 20"



**12AH or 20AH lithium-ion battery**  
Fits in most trunks & through most doorways

**Airline & Cruise Approved**



**MOTOR WARRANTY**

**IQ-8000**  
17.5in seat space



**FRAME WARRANTY**



**MOTOR WARRANTY**

**IQ-8000 PLUS**  
20in seat space



**FRAME WARRANTY**

