

Jazzy[®] carbon HD

Weighs just
32lbs

400lbs
weight capacity



FDA Class II
Medical Device¹

Folds & unfolds
for easy transport



Quick disconnect joystick
with LCD display



Full front &
rear suspension



FDA Class II
Medical Device¹



FDA Class II
Medical Device¹

Heavy-Duty Travel Mobility!

The Jazzy[®] Carbon HD combines a 400-lb. weight capacity with a lightweight mobility option. This power wheelchair weighs just 32.2 lbs. without the battery but doesn't compromise on quality, safety, and the consistent performance you've come to expect from Jazzy[®]. It's easy to fold, easy to pick up, and easy to store and transport. Experience the very best in heavy-duty travel mobility.



FDA Class II
Medical Device*

HCPSC Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- 2 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Includes standard seat weight.

* Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

** Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Consult your airline prior to flight to ensure your lithium battery is permitted.

HCPSC codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

- Unit weighs 32.2 lbs. without the battery
- 400-lb. weight capacity
- Durable carbon fiber frame
- Front and rear suspension
- Top speed of up to 3.7 mph
- Folded dimensions: 31.3" L x 26.4" W x 13 H
- Easy fold design that locks into folded position
- Flip up lightweight carbon fiber arms
- Padded seat back with removable cushion
- Extra-wide adjustable foot platform 17.25" W x 5.9" D
- 20 Ah lithium battery** charges quickly (Standard)
- Storage pocket on back of seat
- Quick disconnect LCD joystick with digital display, USB-C port & puddle light
- Mesh cup holder, cushion and adjustable lap belt, (standard)

SPECIFICATIONS

WEIGHT CAPACITY	400 lbs.	
MAXIMUM SPEED ¹	Up to 3.7 mph	
SEAT DIMENSIONS ³	Width: 20"	Length: 20"
TURNING RADIUS ³	36.4"	
GROUND CLEARANCE ³	1" at motor	
BATTERY WEIGHT ⁴	7.5 lbs	
RANGE PER CHARGE ^{1,2}	Up to 12.4 miles	

COLOR OPTIONS

BLACK CHERRY

NIGHT SKY

OPTIONAL ACCESSORIES

- Additional 12Ah airline-compliant battery**
- Additional 20Ah lithium battery (Coming Soon)
- Battery Travel Bag
- Headrest with Attendant Control Bracket
- Jazzy® Travel Kit (Coming Soon)
- Military Patches
- Simplicity Cushion 20"x 20"
- Elevating Leg Rests (ELRs) (Coming Soon)



Optional Headrest with
Attendant Control Bracket



Optional Jazzy® Travel Kit
(Coming Soon)

WARRANTY

- FRAME: 5-year
- ELECTRONICS: 1-year limited
- DRIVETRAIN: 1-year limited
- BATTERY: 1-year limited