

# Tote Model 003



*WheelChair Carrier*

## Standard Features

- For manual wheelchairs up to 15-1/2" wide when closed
- Lifting the wheelchair onto the carrier is required
- Stationary design does NOT tilt for loading or fold for storage

## Product Specifications

- All steel, & powder coat finish
- 75 lb weight capacity
- Weighs 16 lbs
- Attaches to any ball mount

*Made in the U.S.A*

WheelChair Carrier

Toll Free: (800) 541-3213

E-mail: [customerservice@wheelchaircarrier.com](mailto:customerservice@wheelchaircarrier.com)

[www.WheelChairCarrier.com](http://www.WheelChairCarrier.com)

*Low cost wheelchair carrier!*



Hold down arm  
secures the chair to  
the carrier

