

Mini Electric Lift

Model 117



WheelChair Carrier

Standard Features

- Automatically folds up to stowaway position
- Switch box allows for easy operation of the lift
- Includes 2 retractable ratchets for securing the chair
- Includes wire harness with 30 amp resettable circuit breaker
- **SWING AWAY & BOLT ON RATCHET OPTION AVAILABLE**

Product Specifications

- All steel, & powder coat finish
- 150 lb weight capacity
- 37" long X 22" wide usable platform area
- Weighs 88 lbs
- Can accommodate Class II (1-1/4") or Class III (2") hitch


Made in the U.S.A!

WheelChair Carrier

Toll Free: (800) 541-3213

E-mail: customerservice@wheelchaircarrier.com

www.WheelChairCarrier.com



Transport your lightweight mobility device with ease! This electric lift makes it simple to travel.

