

PAISEEC



One Chair. Three Functions. Unlimited Possibilities.

- 3-in-1 Design: Electric wheelchair, rollator, electric transport chair
- Travel-Ready: Lightweight, foldable, Airline-approved battery
- Smooth & Safe Control: Easy handling for users and caregivers



Standard Colors



Red




Gray

Paiseec W3 - 3-IN-1 Electric Wheelchair

Sort	Features	W3
Vehicle Performance	Maximum Speed	3.7 mph
	Range Per Charge	16 miles
	Maximum Climbing Angle	8°
	Maximum Obstacle crossing height (stepped)	1.57"
	Weight Capacity	300 lbs
Battery Specifications	Battery Specifications	24V 12Ah lithium battery
	Battery Weight	8.2 lbs
	Battery Detachable	Yes
Vehicle Specifications	Ground to Top of Seat	19.3"
	Ground to top of back of seat	33.3"
	Weight (excluding battery)	49.6 lbs
	Front Wheels	8-inch universal wheel
	Rear Wheels	9-inch pneumatic tires
	Ground Clearance	1.57"
Seat Specifications	Seat Dimension (W*D)	19.6"*16.7"
Electronic System	Drive System	"PAI" Smart Safety Riding System
	Motor Specifications	180W*2 brushless motor
	Brake System	Dual brake system: Electromagnetic brake & electronic brake
	Gyroscope	Yes
Storage	Basket	Supported



 TSA Approval



Optional Accessories



Dustproof Cover



Durable Tote Bag



Insulated Cup Holder



Universal Cane Holder