






### Standard Fabrics:




Calypso (TCP)



Camel (TCM)



Hazelnut (THA)



Shiraz (TSZ)



Sterling (TST)

This fabric has a modern, updated look with rich color options. It features our exclusive SofTouch™ with Stain Defense coating.

\* An upcharge will apply for all crossover fabrics.

## EZ Sleeper Slim PR763-MED



*Comfort Zone*<sup>™</sup>  
by GOLDEN

Experience comfort like never before with the EZ Sleeper Slim by Golden. This model offers a modern style track arm design, while the plush back, chaise pad and supportive pocketed coil springs provide hours of relaxing comfort. The EZ Sleeper Slim comes standard with our 2-Zone Comfort Zone positioning system, which features our patented MaxiComfort® positioning operated by our programmable AutoDrive™ hand control.

### Standard Features:

- Patented MaxiComfort® Positioning
- Programmable AutoDrive™ Hand Control
- Extra Pocket
- SmartTek® Diagnostic System



# MAXICOMFORT® Series PR763

Model Number	PR763-MED
Size	Medium
Class II Medical Device	Yes
Weight Capacity	375 LB
Positioning	MaxiComfort
Comfort Zone(s)	2
Twilight Technology	No
ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	No
Arm Rest Covers	No
Overall Width	33"
Overall Height	45"
Seat Width	23"
Floor to Top of Seat	20"
Seat Depth	20.5"
Seat to Top of Back	30"
Seat to Top of Arm	8"
Distance Required from Wall Reclined	23"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 YR
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 YR
Warranty on Electrical Labor	1 YR
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
Heat Wave Technology Option	No
Articulating Headrest Option	No
Power Adjustable Lumbar Option	No
Foot Rest Extension Option	Yes
Straight Lift Option	Standard
USB Charging Port	No
Battery Backup	9 volt
SmartTek® Diagnostic System	Standard
Place of Origin	USA



Shown in Brisa® Cream\*

All Golden lift chairs are Class II Medical Devices.




## Comfort Zone™ by GOLDEN

**Comfort Zones 1 & 2** work together to create individualized wellness positioning and comfort. Each Comfort Zone can be adjusted independently by using the exclusive AutoDrive hand control.



**WE PATENTED COMFORT™**



Scan with your smartphone camera to learn more about Comfort Zone by Golden.



Phone: 800-624-6374  
Fax: 800-628-5165  
www.goldentech.com



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with  
UL 73 and  
CSA C22.2 No. 68



www.facebook.com/goldentechnologies

SHT-PR763-010824