






Standard Fabrics:




Calypso (TCP)



Camel (TCM)



Hazelnut (THA)



Shiraz (TSZ)



Sterling (TST)

This fabric has a modern, updated look with rich color options. It features our exclusive SofTouch™ with Stain Defense coating.

* An upcharge will apply for all crossover fabrics.

EZ Sleeper Slim PR764



Comfort Zone™
by GOLDEN

The EZ Sleeper Slim lift recliner features a modern look with a tailored back and stylish track arm design, while the plush chaise pad and pocketed coil spring seat provide hours of unmatched comfort. Positioning options are endless with our patented 5-Zone Twilight Tilt Technology and MaxiComfort® positioning, adjustable headrest and lumbar support – all controlled effortlessly by our exclusive programmable AutoDrive™ hand control.



FDA Class II Medical Devices designed by Golden Technologies are intended to assist persons with chronic physical disabilities.

Standard Features:

- Patented Twilight Tilt Technology
- Patented MaxiComfort® Positioning
- Adjustable Headrest & Lumbar
- Programmable AutoDrive™ Hand Control
- Extra Pocket
- SmartTek® Diagnostic System



MAXICOMFORT® Series PR764

Model Number	PR764
Size	Medium
Class II Medical Device	Yes
Weight Capacity	400 LB
Positioning	MaxiComfort
Comfort Zone(s)	5
Twilight Technology	Yes
ZG+	Yes
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	No
Arm Rest Covers	No
Overall Width	33.5"
Overall Height	45.5"
Seat Width	22.5"
Floor to Top of Seat	20"
Seat Depth	22"
Seat to Top of Back	30"
Seat to Top of Arm	7"
Distance Required from Wall Reclined	21.5"
Extension Reclined from Back of Seat to End of Foot Rest	42"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 YR
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 YR
Warranty on Electrical Labor	1 YR
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
Heat Wave Technology Option	Yes
Nirvana with Heat Wave Option	Yes
Articulating Headrest Option	Standard
Power Adjustable Lumbar Option	Standard
Foot Rest Extension Option	No
Straight Lift Option	Standard
USB Charging Port	Yes
Battery Backup	Lithium-Ion
SmartTek® Diagnostic System	Standard
Place of Origin	USA



Shown in Redwood

Comfort Zone™ by GOLDEN



Comfort Zones 1 - 5 This positioning system incorporates all five adjustable comfort zones to create the ultimate wellness comfort experience with our patented Twilight Technology. Each Comfort Zone can be adjusted independently by using the exclusive AutoDrive™ hand control.



WE PATENTED COMFORT™



Scan with your smartphone camera to learn more about Comfort Zone by Golden.



Phone: 800-624-6374
Fax: 800-628-5165
www.goldentech.com



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with
UL 73 and
CSA C22.2 No. 68



www.facebook.com/goldentechnologies

LIT-SHT-PR764, Rev. A 060425