

strive ez

ADAPTIVE STROLLER

Everyday Comfort. Limitless Possibilities.

The Strive EZ adaptive stroller combines comfort, support, and ease of use—so every outing is a better one.



EMPOWERING CHILDREN TO MOVE FORWARD

The Strive EZ is a lightweight and easy-to-fold adaptive stroller designed with travel convenience in mind. It provides upper body support, posture maintenance, and comfortable full-body positioning.

It is designed to meet the diverse needs of children with disabilities and offers a wide range of positioning features and optional accessories. Available in four sizes, it accommodates growth with adjustable seat depth. The full-width tilting push bar helps caregivers find a more comfortable pushing angle.



Folds

Easy storage and transport with secure tie-closure strap.





SPECIFICATIONS

	STZ1200	STZ1400	STZ1600	STZ1800
Back Height:	24"	24"	24"	24"
Seat Depth:	9" - 13"	10" - 14"	11" - 15"	12" - 16"
Seat-to-Back Angle:	80° - 95°	80° - 95°	80° - 95°	80° - 95°
Seat Angle:	20°	20°	20°	20°
Seat Width:	12"	14"	16"	18"
Seat-to-Footplate Height:	4" - 20"	4" - 20"	4" - 20"	4" - 20"
Seat-to-Floor Height at 90° (rear/front):	19"/22"	19"/22"	19"/22"	19"/22"
Overall Width (outside rear wheels):	20.5"	22.5"	24.5"	26.5"
Product Weight:	33 lb	34 lb	35 lb	36 lb
Weight Capacity:	75 lb	100 lb	170 lb	250 lb
Weight Capacity for Transit:	75 lb	100 lb	165 lb	165 lb
Front Wheels:	8"	8"	8"	8"
Rear Wheels:	10"	10"	10"	10"
Color:	Red	Red	Blue	Blue

ACCESSORIES

01

Headrest Extension

ST8112 ST8114
ST8116 ST8118

Extends back height and support for taller individuals.



02

Padded Winged Headrest

ST8801

Adjustable padded headrest provides comfortable midline positioning, promoting optimal posture and alignment.



03

Canopy

ST8100S ST8100L

Extendable canopy with a clear flap-covered window. Provides protection from the rain and sunlight.



04

Bus Transit Tie-Downs

ST8802

Secures stroller for safe transportation in a bus/van.



05

Anti-Tippers

ST8804

Prevents backward tipping and enhances stability and security.



06

Upper Extremity Support

ST8805

Depth adjustable, multi-functional cle tray is easy to clean and fits all sizes.



07

Basket

ST8803

Underseat storage basket.

