



smirthwaite

Together we
transform lives

Manufacturer Specification Book

Changing & Showering Product Line

(USA) Version: v1.0 • Date: 2025-12-12

Source: USA Changing & Showering Catalog



Contents



1. Easi-Lift Changing Bench
2. Hi-Riser (Dry & Bath Rooms)
3. Hi-Riser Shower Stretcher
4. Mobi Changer
5. Mobi Heavy Duty Changer
6. Mobi Shower Trolley
7. Shower Stretcher (Fixed Height)
8. Easi-Lift Shower Stretcher
9. Warranty & Contact



Easi-Lift Changing Bench

Product Overview

Key Features:

- No floor fixing required
- Electrically operated
- Folds away to save space
- Stowaway safety rail
- Weight capacity ~330 lbs
- Travel range ~17.7–37.4 in height



Measurement	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6	Size 7	Size 8
Code	3081	3082	3083	3084	3085	3086	3087	3088
Table dimension length (inches)	47.2	51.8	55.1	59.0	62.9	66.9	78.7	82.6
Table dimension width (inches)	21.6	21.6	23.6	23.6	25.5	25.5	25.5	25.5
Table dimension width including safety rail (inches)	34.4	34.4	36.4	36.4	38.3	38.4	38.4	38.4
Travel height range (inches)	17.7 - 37.4	17.7 - 37.4	17.7 - 37.4	17.7 - 37.4	17.7 - 37.4	17.7 - 37.4	17.7 - 37.4	17.7 - 37.4
Wall frame width (inches)	27.6	27.6	27.6	27.6	27.6	27.6	27.6	27.6
Wall frame depth (inches)	4.3	4.3	4.3	4.3	4.3	4.3	4.3	4.3
Wall frame height (inches)	63.7	63.7	63.7	63.7	63.7	63.7	63.7	63.7
Safety side height (inches)	6.5	6.5	6.5	6.5	6.5	6.5	6.5	6.5
Stored width (inches)	12.0	12.0	12.0	12.0	12.0	12.0	12.0	12.0
Max user weight (lbs)	330	330	330	330	330	330	330	330
Product weight (lbs)	176	180	185	189	194	198	207	211


smirthwaite

Hi-Riser (Dry & Bath Rooms)

Product Overview

Key Features:

- Electrically operated
- Wall and floor mounted
- Weight capacity ~441 lbs
- Travel range ~24 in height
- Frames: cream powder coated (dry rooms) or stainless steel (bath rooms)
- Folds away to save space
- Easy to clean padded surfaces
- Stowaway side rail



Measurement	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6	Size 7	Size 8	Size 9
Code for mild steel powder coated	3041	3042	3043	3044	3045	3046	3047	3048	3049
Code for stainless steel	3051	3052	3053	3054	3055	3056	3057	3058	3059
Table dimension length (inches)	47.2	51.1	55.1	59	62.9	66.9	70.8	78.7	82.6
Table dimension width (inches)	21.6	21.6	23.6	23.6	25.6	25.6	25.6	25.5	25.6
Table dimension with with safety guard (inches)	34.4	34.4	36.4	38.4	38.4	38.4	38.4	38.4	38.4
Travel range (inches)	15.7 - 39.4	15.7 - 39.4	15.7 - 39.4	15.7 - 39.4	15.7 - 39.4	15.7 - 39.4	15.7 - 39.4	15.7 - 39.4	15.7 - 39.4
Wall frame width (inches)	27.6	27.6	27.6	27.6	27.6	27.6	27.6	27.6	27.6
Wall frame depth (inches)	4.3	4.3	4.3	4.3	4.3	4.3	4.3	4.3	4.3
Wall frame height (inches)	61.0	61.0	61.0	61.0	61.0	61.0	61.0	61.0	61.0
Safety guard height (inches)	6.5	6.5	6.5	6.5	6.5	6.5	6.5	6.5	6.5
Stored width (inches)	12	12	12	12	12	12	12	12	12
Max user weight (lbs)	440	440	440	440	440	440	440	440	440
Product weight (lbs)	176	180	185	189	194	198	202	207	212

smirthwaite

Hi-Riser Shower Stretcher

Product Overview

Key Features

- Electrically operated handset (safe in wet rooms)
- Two-in-one changing & showering
- Drainage pipe & fasteners included
- Wall and floor mounted
- Stainless steel construction
- Weight capacity ~331 lbs
- Travel range ~22.4 in height
- Folds away; only ~11.8 in from wall when stored or mounted



Measurement	Size 1	Size 2
Code for stainless steel	3063	3067
Table dimension length (in)	55.1	70.8
Table dimension width (in)	25.5	25.5
Table dimension width with safety guard (in)	43.3	43.3
Travel range (in)	15.7 - 39.4	15.7 - 39.4
Wall frame width (in)	27.6	27.6
Wall frame depth (in)	4.3	4.3
Wall frame height (in)	61	61
Safety guard height (in)	6.5	6.5
Stored width (in)	12	12
Max user weight (lbs)	331	331
Product weight (lbs)	176	187



smirthwaite

Mobi Changer

Product Overview

Key Features

- Mobile; no installation required
- Frames: cream powder coated (dry rooms) or stainless steel (wet rooms)
- Available lengths: 55.1 in or 70.9 in
- Weight capacity ~500.4 lbs
- Travel range ~37.4 in height
- Optional angled back
- Unique stowaway side rails



Measurement	Size 1	Size 2	Size 3	Size 4
Code for mild steel powder coated	3121	3122	3123	3124
Code for stainless steel	3141	3142	3143	3144
Table dimensions length (inches)	55.1	62.9	70.9	78.7
Table dimensions width (inches)	29.5	29.5	29.5	29.5
Travel range (inches)	18.5 - 37.4	18.5 - 37.4	18.5 - 37.4	18.5 - 37.4
Safety guards width (inches)	33.1	33.1	33.1	33.1
Safety guards height (inches)	6.5	6.5	6.5	6.5
Safety guards operational space (inches)	7.9	7.9	7.9	7.9
Castor centres length (inches)	42.9	42.9	42.9	42.9
Castor centres width (inches)	27.6	27.6	27.6	27.6
Nominal clearance (inches)	7.5	7.5	7.5	7.5
Max user weight (lbs)	500.4	500.4	500.4	500.4
Product weight (lbs)	189	194	198	220

Mobi Heavy Duty Changer

Product Overview

Key Features

- Mobile; no installation required
- Electronically operated
- Dry rooms (cream powder-coated frame)
- Length 78.7 in
- Weight capacity ~1000 lbs
- Travel range ~36.2 in height
- Unique stowaway side rails
- Wall-mounted charger & 2 batteries included
- Rooms (cream powder-coated frame)



Measurement	Size 1
Code for mild steel powder coated	3161
Table dimensions length (Inches)	78.7
Table dimensions width (Inches)	37.8
Travel range (Inches)	21.7 - 36.2
Safety guards width (Inches)	41.3
Safety guards height (Inches)	6.5
Safety guards operational space (Inches)	7.9
Castor centres length (Inches)	59.1
Castor centres width (Inches)	31.5
Nominal clearance (Inches)	5.9
Max user weight (lbs)	1000
Product weight (lbs)	220



Mobi Shower Trolley

Product Overview

Key Features

- Mobile; no installation required
- Stainless steel frame
- Drainage channel & pipe
- Soft mattress & wedge pillow
- Length 70.9 in
- Weight capacity ~485 lbs
- Travel range ~18.3–36.4 in height
- Unique stowaway side rails



Measurement	Size 1
Code for stainless steel	3151
Table dimensions length (inches)	70.9
Table dimensions width (inches)	29.5
Travel range (inches)	18.3 - 36.4
Safety guards width (inches)	33.1
Safety guards height (inches)	7.7
Safety guards operational space (inches)	7.9
Castor centres length (inches)	42.9
Castor centres width (inches)	27.6
Nominal clearance (inches)	7.5
Max user weight (lbs)	500.4
Product weight (lbs)	189

smirthwaite

Shower Stretcher (Fixed Height)

Product Overview

Key Features

- No floor attachments required
- Ideal for wet rooms
- Stainless steel frame with powder-coated finish
- Permanently fixed at desired height
- Available lengths: 55.1 in & 74.8 in
- Weight capacity ~440.9 lbs
- Folds away; wipe-clean surfaces
- Includes safety guard
- Optional integral vinyl headrest



Measurement	Size 1	Size 6
Code stainless steel powder coated w/out headrest	3171	3176
Code stainless steel powder coated with headrest	3271	3276
Table dimension length (inches)	55.1	74.8
Table dimension width (inches)	25.6	25.6
Wall frame width (inches)	26.2	26.2
Wall frame depth (inches)	3.1	3.1
Stored width (inches)	9.4	9.4
Max user weight (lbs)	440	440
Product weight (lbs)	165	176

Easi-Lift Shower Stretcher

Product Overview

Key Features

- No floor attachments required
- Ideal for wet rooms
- Stainless steel frame
- Travel range 17.7 - 37.4 inches
- Available lengths: 55.1 in & 74.8 in
- Weight capacity of 441 lbs
- Folds away; wipe-clean surfaces
- Includes safety guard
- Optional integral headrest



Measurement	Size 6	Size 1
Stainless steel without headrest	3196	3191
Table dimensions length (in)	74.8	55.1
Table dimensions width (in)	25.6	25.6
Travel height range (in) - dependent on installed height	17.7 - 37.4	17.7 - 37.4
Wall frame width (in)	29.5	29.5
Wall frame depth (in)	4.3	4.3
Wall frame height (in)	51.9	51.9
Safety guard height (in) - optional extra	5.7	5.7
Stored width (mm)	10.6	10.6
Max user weight (lbs)	441	441
Product weight (lbs)	192	180

Warranty

Standard Warranty:

2 years on all products. For details or questions contact your local representative.

Together we
transform lives

Contact

212 541 2455

info@smirthwaite.com

smirthwaite.com

