TRANSITIONS®

Threshold ramps, mats, and plates



EZ-ACCESS[®] knows that a gentle slope can make a world of difference. TRANSITIONS[®] threshold solutions allow for a smooth ground-to-sill transition over doorway thresholds throughout your home. Reduce preventable trip hazards and make your home more comfortable and safe for mobility users.



The TRANSITIONS® Advantage

- Sturdy, seam-free surface makes TRANSITIONS ramps, mats, and plates the safe solution for wheelchair, scooter, walker, and cane users
- Attractive designs are available in multiple styles and materials to suit your needs and taste
- Slip-resistant tread and surfaces
- Durably designed for years of continued use
- Quick and easy setup and removal installs in minutes so there's little-to-no disruption for you and your family
- Ideal for doorways, sliding glass doors, and raised landings
- Portable, convenient, and durable
- The entire line is backed by our lifetime warranty



Visit ezaccess.com for product specifications

TRANSITIONS® Modular Entry Ramp

- Modular design with interchangeable extensions can accommodate thresholds from 3/4" to 6-1/4" high
- Installation is easy with our hidden hinge ... simply add and remove extensions to accommodate your needs
- Pivoting flap attaches to the top of the ramp to provide 1/4" of additional adjustability
- Made of aluminum and features an extruded, slip-resistant surface for superior traction
- 34" wide
- 700 lb. weight capacity



TRANSITIONS[®] Angled Entry Ramp



- Swivel feet accommodate varying heights, angles, and surface types
- Three lengths are available (12", 24", and 36") to accommodate thresholds ranging from 1-3/8" to 5-7/8" high
- Manufactured using high-strength aluminum that will not rust or corrode
- Features an extruded, slip-resistant surface for superior traction
- 36" wide
- 700 lb. weight capacity

TRANSITIONS® Angled Entry Mat

- This "welcome" mat made using recycled rubber is eco-friendly and can be used safely on nearly any surface
- · Stands up well in harsh outdoor conditions
- Features beveled edges for cross-traffic access
- Available in 1-1/2" and 2-1/2" heights and can easily be trimmed or notched if needed
- Offered with or without the recessed "Welcome"
- Black or grey color options
- Usable width of 36" (40"-43" overall)
- 850 lb. weight capacity

customerservice@ezaccess.com

TRANSITIONS® Angled Entry Plate

- Formed aluminum construction for simple, one-step threshold navigation
- Three lengths are available (8", 10", and 12") to accommodate thresholds up to 2" high
- Features an applied, slip-resistant surface for greater traction, even when wet
- 32" wide
- 700 lb. weight capacity



1.800.451.1903 | www.ezaccess.com

Text and images © Homecare Products, Inc. All rights reserved.