















The Forcemech Traveler is a lightweight, foldable travel scooter designed for easy portability and reliable performance. Weighing just 38 lbs, it folds swiftly and can be transported effortlessly, making it ideal for those with active lifestyles. Its advanced LED technology includes an LED dash displaying battery power and speed, along with automatic LED lights and turn signals to enhance visibility and safety. Equipped with a quick-swap battery system for uninterrupted journeys, the Traveler combines sleek design, practicality, and safety, making it the perfect travel companion.

## **Product Features**

- ·Simple and Elegant Design
- Smooth Brushless Motor
- Automatic LED Headlight
- •LED Daytime Running Light
- •LED Rear Signal Flow Light
- •Front Dash LED Info Screen
- Hot-Swap Lithium Battery
- •Foldable in Seconds for Travel
- •38Lbs\* Scooter Weight

WWW.FORCEMECH.COM

1-877-903-6723



Model	Forcemech Traveler / FMTL1G1
Material	Aluminum Alloy
Dimension	40.16 x 18.11 x 35.04
Folded Dimension	40.16 x 18.11 x 16.54
N.W.(w/o battery)	38.36 lbs
Front Wheel Size	7.87 x 1.97
Rear Wheel Size	7.87 x 1.97
Loading Capacity	275 lbs
Motor Rated Power	120W
Motor Output Power (Max)	240W
Battery	24V/10AH (Lithium)
Brake Type	Electromagnetic brake 24V4N.m
Controller	Brushless
Charger	AC100-240V50/60HZ Output: 29.2V-2A/3A
Output	AC100-240V50/60HZ Output: 29.2V-2A/3A
Max Speed	4.35 mph
Slope Holding Performance	9°
Overrun Height	1.57 inches
Trench Crossing Wide	3.94 inches

Climbing Capability	≤9°
Mini Turning Radius	≤47.24 inches
Standard Driving Range	6.21 miles (10Ah); 12.43 miles (20Ah)
Upper Charging Port	Available
Ground Clearance	2.95 inches
Model	Forcemech Traveler / FMTL1G1
Material	Aluminum Alloy
Dimension	40.16 x 18.11 x 35.04
Folded Dimension	40.16 x 18.11 x 16.54
N.W.(w/o battery)	38.36 lbs
Front Wheel Size	7.87 x 1.97
Rear Wheel Size	7.87 x 1.97
Loading Capacity	275 lbs
Motor Rated Power	120W
Motor Output Power (Max)	240W
Battery	24V/10AH (Lithium)
Brake Type	Electromagnetic brake 24V4N.m
Controller	Brushless

WWW.FORCEMECH.COM