

merits

# VISION SPORT P326E

MWD Powerbase

K0822 / K0823



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

# Features

- True mid-wheel: 20" turning Radius
- Robust independent suspension can absorb up to 250 lbs. on the footplate without tipping.
- Premium, semi-reclining ergonomic seat, up to 135°
- 10" Solid drive tires for superior traction
- Dynamic LiNX 40A controller for convenient Bluetooth programming and diagnostics
- 16" - 20" customizable seat depth

# Options

- Unoccupied transit tie down kit
- Cup holder
- USB charger
- Oxygen tank holder
- Walker holder
- 16" - 20" seat width option
- Swing-away joystick
- Elevating leg rests (K0195)
- Foldable rear basket
- Cane & crutch holder



# Specifications

Product#	P326
Weight Capacity	300 lbs.
Drivetrain	MWD
Seat Type	Captain or Pan
Seat Width	16" - 20"
Seat Depth	16" - 20"
Backrest Height	18" - 20"
Seat-to-Floor Height	17" - 21"
Overall Length	34"
Overall Width	24"
Overall Height	48.5"
Turning Radius	20"
Ground Clearance	2"
Total Weight	184 lbs.
Wt of Heaviest Pc.	90 lbs.
Max Speed	5 mph
Range Per Charge	up to 18 mi.
Battery Size (2 required)	33Ah, U1 AGM
Motor	DC24V / 160W
Controller	Dynamic LiNX 40A
Charger	4A off-board
Gradient	6°
Front Wheels	6" - solid
Drive Wheels	10" - solid
Rear Wheels	6" - solid
HCPCode	K0822 / K0823
Color Options	<div><div></div> Red</div> <div><div></div> Black</div>

Features MK batteries



**merits** *Enhancing life through mobility.*

©Merits Health Products, Inc.  
4245 Evans Avenue Fort Myers, FL 33901  
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com

