

Nitro® Sprint Rollator with Accessories

HCPCS:



The sleek Nitro Sprint Rollator is designed to provide comfort, convenience, and control at every turn. The rollator sports a unique braking system that reduces tire wear and has an integrated slow-down brake for safer maneuverability over diverse terrain. When not in use, the rollator can be folded side to side with one hand and locked in place for easy transport. Other features include multiple storage bag mounting points and a wider backrest with a window for improved comfort and visibility.

- Built-in accessory mount easily attaches a cup, cane or smartphone holder. *ALL 3 ACCESSORIES INCLUDED! (Cup Holder, Phone Holder, Cane Holder)
- Contoured backstrap with viewing window for improved comfort and visibility.
- The patented frame design provides added strength and durability with a 350lb weight capacity.
- Integrated slow-down brake controls speed when navigating diverse terrain.
- The enhanced braking system reduces tire wear and allows for easy replacement of brake pads and cables.
- Improved tire tread for better stopping ability.
- Built-in brake cable for improved durability and safety
- Large 10" front wheels provide optimal stability and rolling comfort.
- Folds side to side with one hand and remains locked for easy transport and storage.
- A padded seat provides added comfort while sitting.
- Storage bag mounts in multiple locations for easy access.

PRODUCTS	Description
102662BK-ACC	Nitro Sprint Rollatr Bk+Cup+Cell+CanePKG
102662BK-H-ACC	Nitro Sprint Rollatr BkH+Cup+Cell+CnePKG
102662BK-T-ACC	Nitro Sprint Rollatr BkT+Cup+Cell+CnePKG
102662BL-ACC	Nitro Sprint Rollatr Bl+Cup+Cell+CanePKG
102662BL-H-ACC	Nitro Sprint Rollatr BlH+Cup+Cell+CnePKG
102662BL-T-ACC	Nitro Sprint Rollatr BlT+Cup+Cell+CnePKG
102662RD-ACC	Nitro Sprint Rollatr Rd+Cup+Cell+CanePKG
102662RD-H-ACC	Nitro Sprint Rollatr RdH+Cup+Cell+CnePKG
102662RD-T-ACC	Nitro Sprint Rollatr RdT+Cup+Cell+CnePKG

SPECIFICATIONS	Description
Warranty	Frame: Limited Lifetime, Brake Handle & Housing: 5 Years
Product Weight	Hemi/Standard: 19 lb Tall: 19.5 lb
Handle Height	Hemi: 29-33 Standard: 30-36.5 Tall: 32-39.4
Weight Capacity	350 lbs
FireRating	Back and Seat fabric TB117-2013 section 1
Seat To Floor Height	Hemi: 18" Standard: 21" Tall: 23"

SPECIFICATIONS**Description**

Dimensions Assembled	27.5" (L) X 24" (W) X 30 -36.5" (H)
Width inside back legs	19"
CA Prop 65	no