

Jazzy[®]

POWER CHAIRS



America's Favorite Power Chair



The Advantages of the Jazzy® Power Chair



When choosing a Jazzy Power Chair, specific requirements and conditions can affect your process...

The Jazzy® Power Chair is consumer-inspired and engineered to provide you with consistent high-performance operation you can count on. With a full line of high-performing Jazzy Power Chair models available and many accessories and options to choose from, there's a Jazzy to meet your unique needs and lifestyle.

Indoor Maneuverability

Compact and mid-wheel design offer Jazzy's legendary maneuverability

Increased Weight Capacity

Most full-size power chairs offer increased weight capacities

Stability

Active-Trac® ATX Suspension for enhanced performance and comfort over more varied terrain

Easy Portability

Portable Jazzy models offer easy disassembly for convenient transport

Indoor/Outdoor Versatility

Jazzy power chairs are great for your lifestyle, indoors and out

Speed and Performance

Mid-wheel drive, Active-Trac ATX Suspension for increased speed and performance

What's Inside...

Models:

Non-Coded.....4-6

PDAC-Coded....7-11

Accessories.....12

Innovative Technology.....13

Specifications.....14-15

Jazzy Air®

- 300 lb. weight capacity
- Speeds up to 4 mph non-elevated and 3.5 mph elevated with Extra Stability Technology
- Up to 15 miles with U-1 batteries and 18.4 miles w/40 Ah batteries per charge range
- Mid-wheel 6 technology with Active-Trac® ATX Suspension
- Choice of 16" x 16"-18", 18" x 18"-20" and 20" x 20"-22" depth-adjustable high-back reclining sport seat



Recipient of the 2016
HME Innovative Retail
Product Award

Available Colors:



Go Chair®

- Ultimate in portability
- 300 lb. weight capacity
- Faster speed up to 3.7 mph
- Up to 8.7 miles per charge range
- Dual swing-away, under seat storage bins
- One hand feather-touch disassembly



Go Chair®
shown disassembled

Available Colors:



Front arm mount



Under arm mount

Air Button

The **Air Button** allows you to experience 10" of power adjustable seat height while on the move, and can be mounted in multiple locations for easier access.



Storage

Dual, swing-away under seat compartments for convenient storage of personal items.



Transport

Feather-touch disassembly for easy, lightweight travel options.



Jazzy® Sport 2

- 300 lb. weight capacity
- Speeds up to 4 mph
- Up to 13 miles per charge range
- Front-wheel drive for enhanced maneuverability

Available Colors:



Jazzy
Red

Jazzy
Blue



Jazzy® Elite ES-1

- 300 lb. weight capacity
- Speeds up to 4 mph
- Up to 12 miles per charge range
- Front-wheel drive

Available Colors:



Jazzy
Red

Jazzy
Blue



Jazzy® Sport Portable

- 300 lb. weight capacity
- Speeds up to 4 mph
- Up to 11 miles per charge range
- Easily disassembles into 5 pieces for convenient transport

Available Colors:



Jazzy
Red

Jazzy
Blue



Jazzy® Elite ES Portable

- 300 lb. weight capacity
- Speeds up to 4 mph
- Up to 11 miles per charge range
- Easily disassembles into 5 pieces for convenient transport
- Front-wheel drive

Available Colors:



Jazzy
Red

Jazzy
Blue



Jazzy® Elite ES

- 300 lb. weight capacity
- Speeds up to 4 mph
- Up to 15 miles per charge range
- Front-wheel drive

Available Colors:



Jazzy® Select 6

- 300 lb. weight capacity
- Speeds up to 4 mph
- Up to 15 miles per charge range
- Mid-wheel 6 technology with Active-Trac® ATX Suspension

Available Colors:



Jazzy® Select Elite

- 300 lb. weight capacity
- Speeds up to 4 mph
- Up to 15 miles per charge range
- Front-wheel drive

Available Colors:



Jazzy® Elite 14

- 300 lb. weight capacity
- Speeds up to 4 mph
- Up to 9.8 miles per charge range
- 14" drive wheels for enhanced outdoor performance

Available Colors:






Jazzy® 600 ES

- 300 lb. weight capacity
- Speeds up to 4 mph
- Up to 17 miles per charge range
- Mid-wheel 6 technology with Active-Trac® ATX Suspension

Available Colors:



Jazzy® 614HD

- 450 lb. weight capacity
- Speeds up to 4.3 mph
- Up to 11 miles per charge range
- Mid-wheel 6 technology with Active-Trac® ATX Suspension

Available Colors:



Jazzy® Elite HD

- 450 lb. weight capacity
- Speeds up to 4 mph
- Up to 9.8 miles per charge range
- Front-wheel drive

Available Colors:



Jazzy® 1450

- 600 lb. weight capacity
- Speeds up to 4.5 mph
- Up to 12 miles per charge range
- Aggressive front-wheel drive performance

Available Colors:



Accessories

A wide selection of accessories are available to personalize your Jazzy® and make it uniquely yours.



Cup Holder



Jazzy Weather Cover



Saddle Bag



XLR USB Charger



Rear Mount
Cane/Crutch Holder



Walker Holder



Elevating Leg Rests



Cell Phone Holder



Swing Away
Leg Rests



Heel Loops



Safety Flag



Swing-Away
Joystick



Rear Basket



Oxygen Holder



Lap Belt
(additional sizes)



Military
Patches



Innovative Technology

Jazzy Air® features seating that elevates 10" in just 16 seconds to primary standing height, allowing exceptional access to your environment. Extra Stability Technology mechanically stabilizes the power base when elevated, creating ultimate security to safely drive up to 3.5 mph at walking height.

Mid-Wheel 6® technology with Active-Trac® ATX Suspension offers the highest levels of performance and maneuverability while easily and smoothly climbing transitions and obstacles. This patented technology allows the active front casters to both lift when climbing a threshold, while also dropping down for stability when descending.

In-line motor technology offers enhanced torque, battery range and reliability for the most dependable performance.

Consumer-inspired and real-world proven, Jazzy® offers innovative technologies you can count on.

DISCLAIMER: Not all accessories are available for all products. Please speak with your provider for details and availability.

Specifications

	MODEL	Weight Capacity	Turning Radius ¹	Width ¹	Length ^{1,2}	Maximum Speed ³	Ground Clearance ¹	Front Wheels	Drive Wheels	Rear Wheels	Drivetrain	Suspension Type	High-Back Seat (max. size)	High-Back Seat-to-Floor Range	Synergy [®] Seating	Synergy [®] Seat-to-Floor Range	Power Elevating Seat ⁷	Standard Electronics	Battery Size (2 required)	Standard Battery Charger	Per-Charge Range ^{5,6}	Battery Weight ⁹	Weight of Heaviest Piece ⁵	Standard Seat Weight	Standard Seat Size	HCPCS Codes	Warranty
NON-CODED	GO-CHAIR	300 lbs.	25.5"	22.5"	28.5"	Up to 3.7 mph	2" @ frame	6" solid casters	9" solid	3" solid anti-tips	2-motor, RWD	Limited	N/A	N/A	N/A	N/A	N/A	40A, PG GC3	18 AH, 12V deep-cycle	2A, off-board	Up to 8.7 miles	14 lbs. each	36 lbs.	30 lbs. compact, foldable	18" x 17"	K0899	5-YR. LIMITED: Frame 13 MO. LIMITED: electronics/drive motors
	JAZZY AIR [®]	300 lbs.	22.5"	25.5"	38.5"	EL= Up to 4 mph NON-EL= Up to 3.5 mph	2.42" @ foot-plate	6" solid casters	10"	6" solid casters	2-motor, in-line, MWD	Active-Trac [®] (ATX)	20"W x 20" x 22"D	23"-33"	N/A	N/A	Std.	40A, PG VR2	12 V, U-1	3.5A, off-board	Up to 15 miles ⁹	24.5 lbs. each	100 lbs.	41 lbs.	18" x 18" -20"	K0899	5-YR. LIMITED: Frame 13 MO. LIMITED: electronics/drive motors
	JAZZY [®] SPORT 2	300 lbs.	24.75"	22.75"	34.75"	Up to 4 mph	1.5"	3" solid casters	9"	6" solid casters	2-motor, in-line, FWD	Limited	20"W x 20"D	21.5" x 23.5"	N/A	N/A	N/A	40A, PG GC3	12 V, U-1	3A, off-board	Up to 13 miles	24.5 lbs. each	69 lbs.	43 lbs. comfort, high-bk.	18" x 18" -20"	K0899	5-YR. LIMITED: Frame 13 MO. LIMITED: electronics/drive motors
	JAZZY [®] SPORT PORTABLE	300 lbs.	25.5"	22.75"	34.75"	Up to 4 mph	1.5"	3" solid anti-tips	9" solid	6" solid casters	2-motor, in-line, FWD	Limited	20"W x 20"D	21.5" x 23.5"	N/A	N/A	N/A	40A, PG GC3	12 V, U-1	3A, off-board	Up to 11 miles	24.5 lbs. each	48 lbs.	43 lbs. comfort, high-bk.	18" x 18" -20"	K0899	5-YR. LIMITED: Frame 13 MO. LIMITED: electronics/drive motors
PDAC CODED	JAZZY [®] ELITE ES-1	300 lbs.	23.5"	22.75"	34.75"	Up to 4 mph	0.75"	5" solid anti-tips	9" solid	6" solid casters	2-motor, in-line, FWD	Limited	20"W x 20"D	21.5" x 23.5"	N/A	N/A	N/A	40A, PG GC3	12 V, U-1	3A, off-board	Up to 12 miles	24.5 lbs. each	69 lbs.	43 lbs. comfort, high-bk.	18" x 18" -20"	Group 1: Std. (K0816)	5-YR. LIMITED: Frame 13 MO. LIMITED: electronics/drive motors
	JAZZY [®] ELITE ES PORTABLE	300 lbs.	25.5"	22.75"	34.75"	Up to 4 mph	1.5"	3" solid anti-tips	9" solid	6" solid casters	2-motor, in-line, FWD	Limited	20"W x 20"D	21.5" x 23.5"	N/A	N/A	N/A	40A, PG GC3	12 V, U-1	3A, off-board	Up to 11 miles	24.5 lbs. each	48 lbs.	43 lbs. comfort, high-bk.	18" x 18" -20"	Group 2: Std., Portable (K0820/K0821)	5-YR. LIMITED: Frame 13 MO. LIMITED: electronics/drive motors
	JAZZY [®] ELITE ES	300 lbs.	24.75"	22.75"	34.75"	Up to 4 mph	1.5"	3" solid anti-tips	9" solid	6" solid casters	2-motor, in-line, FWD	Limited	20"W x 20"D	21.5" x 23.5"	N/A	N/A	N/A	40A, PG GC3	12 V, U-1	3A, off-board	Up to 15 miles	24.5 lbs. each	69 lbs.	43 lbs. comfort, high-bk.	18" x 18" -20"	Group 2: Std., (K0822/K0823)	5-YR. LIMITED: Frame 13 MO. LIMITED: electronics/drive motors
	JAZZY [®] SELECT ELITE	300 lbs.	24.75"	22.75"	34.75"	Up to 4 mph	1.5"	3" solid anti-tips	9" solid	6" solid casters	2-motor, in-line, FWD	Limited	20"W x 20"D	21.5" x 23.5"	20"W x 20"D	17.25"-19.25"	N/A	40A, PG GC3	12 V, U-1	3A, off-board	Up to 15 miles	24.5 lbs. each	69 lbs.	43 lbs. comfort, high-bk.	18" x 18" -20"	Group 2: Std., (K0822/K0823)	5-YR. LIMITED: Frame 13 MO. LIMITED: electronics/drive motors
	JAZZY [®] SELECT 6	300 lbs.	21.75"	23.5"	35.5"	Up to 4 mph	1.75"	6" solid casters	10" solid	6" solid casters	2-motor, in-line, MWD	Active-Trac [®] (ATX)	20"W x 20"D	21.5"-23.25"	20"W x 20"D	16.5"-18.5"	Optional	40A, PG GC3 ⁴	12 V, U-1	3A, off-board	Up to 15 miles	24.5 lbs. each	100 lbs.	43 lbs. comfort, high-bk.	18" x 18" -20"	Group 2: Std., (K0822/K0823)	5-YR. LIMITED: Frame 13 MO. LIMITED: electronics/drive motors
	JAZZY [®] ELITE 14	300 lbs.	24"	24.25"	35.5"	Up to 4 mph	3"	3" solid anti-tips	14"	6" solid casters	2-motor, in-line, FWD	Limited	20"W x 20"D	20"-22"	20"W x 20"D	15.75"-17.75"	N/A	50A, Dynamic Shark	12 V, NF-22	5A, off-board	Up to 9.8 miles	37 lbs. each	99 lbs.	51.5 lbs. high-bk.	20" x 20" -18"	Group 2: Std., (K0822/K0823)	5-YR. LIMITED: Frame 13 MO. LIMITED: electronics/drive motors
	JAZZY [®] 600 ES	300 lbs.	20.5"	24.5"	35.5"	Up to 4 mph	2.75"	6" solid OMNI-casters	14"	6" solid casters	2-motor, in-line, MWD	Active-Trac [®] (ATX)	20"W x 20"D	21"-24"	24"W x 24"D	17.25"-19.25"	N/A	50A, Dynamic Shark	12 V, NF-22	5A, off-board	Up to 17 miles	37 lbs. each	124 lbs.	52 lbs. high-bk.	18" x 18"/20" x 20" -18"	Group 2: Std., (K0822/K0823)	LIFETIME. LIMITED: Frame 13 MO. LIMITED: electronics/drive motors
	JAZZY [®] ELITE HD	450 lbs.	24"	24.25"	35.5"	Up to 4 mph	3"	3" solid anti-tips	14" solid	6" solid casters	2-motor, in-line, FWD	Limited	24"W x 24"D	20" x 22"	24"W x 24"D	15.75"-17.75"	N/A	60A, Dynamic Shark	12 V, NF-22	5A, off-board	Up to 9.8 miles	37 lbs. each	99 lbs.	52 lbs. high-bk.	20" x 18" -20"	Group 2: Heavy Duty (K0824/K0825)	5-YR. LIMITED: Frame 13 MO. LIMITED: electronics/drive motors
	JAZZY [®] 614 HD	450 lbs.	22.75"	25.25"	35.5"	Up to 4.3 mph	2.25"	6" solid casters	14"	6" solid casters	2-motor, in-line, MW6 ⁶	Active-Trac [®] (ATX)	24"W x 24"D	20.5" x 23.5"	24"W x 24"D	17.75"-19.75"	N/A	50A, Dynamic Shark	12 V, NF-22	5A, off-board	Up to 11 miles	37 lbs.	122 lbs.	52 lbs. high-bk.	20" x 18" -20"	Group 2: Heavy Duty (K0824/K0825)	LIFETIME. LIMITED: Frame 13 MO. LIMITED: electronics/drive motors
	JAZZY [®] 1450	600 lbs.	26.5"	29.25"	38.25"	Up to 4.5 mph	1.5"	4" solid anti-tips	14" solid	9" solid casters	2-motor, FWD	Swing Arm	32"W x 28"D	20.5" x 22.5"	32"W x 28"D	17.25"-19.25"	N/A	90A, PG VR2	12 V, Group-24	8A, off-board	Up to 12 miles	52 lbs. each	195 lbs.	52 lbs. high-bk.	24" x 24"	Group 2: Very Heavy Duty (K0826/K0827), Single Power Very Heavy Duty (K0839)	K0826 / K0827 Models-LIFETIME. LIMITED: Frame 13 MO. LIMITED: electronics/drive motors K0839 Model-LIFETIME. LIMITED: Frame 2 YR. LIMITED: electronics 18 MO. LIMITED: drive motors

- Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.
- Without front riggings.
- Varies with user weight, terrain, battery type, battery condition and tire pressure.
- 40A, PG 6-key VR2 required with Power Elevating Seat Option.
- Does not include seat or battery weights.
- Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Testing conducted at maximum weight capacities. This specification can be subject to a variance of +10%, -5%.
- Maximum weight capacity is reduced to 400 lbs. (181 kg) on 18" deep seat configurations.
- Battery weight may vary based on manufacturer.
- Up to 18.4 miles with 40 Ah batteries

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.



Authorized Pride® Provider

D E A L E R S T A M P

(US) 800-800-8586 • (CAN) 888-570-1113 • pridemobility.com •    