

PRIDE<sup>®</sup>  
POWER  
RECLINER  
COLLECTION



“

I LOVE MY PRIDE LIFT CHAIR!  
I CAN FULLY RECLINE AND BE  
COMPLETELY COMFORTABLE!  
AND, THE FABRIC MATCHES  
MY LIVING ROOM PERFECTLY!

– GEORGE R.



SCOTT MEUSER  
CHAIRMAN & CEO  
PRIDE MOBILITY®

With over three decades in the power lift recliner industry, we know that our business is about you, our customer. We're humbled when you choose to bring one of our power lift recliners into your home, knowing that it must fit your lifestyle and decor.

Pride® Power Lift Recliners is the culmination of your direct feedback. You wished for luxurious comfort, exceptional quality, state-of-the-art features and elegant designs. Over a 30-year period, we listened and refined, then listened some more. The result is the Pride line of power lift recliners designed for you, your home and your lifestyle.

We understand that bringing a power lift recliner into your home isn't just a purchase, but an investment – toward your comfort, home décor, and lifestyle. We deeply value that you, our customer, have enabled us to meet those vital goals through Pride Power Lift Recliners.

A handwritten signature in black ink that reads "Scott Meuser". The signature is fluid and cursive, with a long horizontal stroke at the end.



RECLINE POSITION 1.



2. UPWARD LIFT



FOOTREST EXTENSION 4.



USB REMOTE 5.



BACKREST 6.

- 5. The USB remote adjusts all functions independently to help you find the perfect position.
- 6. The backrest moves up and down independently to help you find the perfect position.



### FUNCTIONS AND FEATURES

- 1. The full recline position is ideal for napping or relaxing.
- 2. The upward lift offers convenient stand assist for getting in and out of the recliner.
- 3. Dual side pockets are ideal for storing remotes, books or magazines.
- 4. The footrest extension allows you to relax and stretch out completely.

## ELEGANT COMFORT AND RELAXATION



\*QUAD SIDE POCKETS



7. INFINITE POSITION

### STANDARD FEATURES

- 7. Infinite position recline enables the backrest and footrest to move independently from one another so that a variety of positions can be reached.

\*Only available on the Oasis Collection.




MY PRIDE RECLINER GIVES ME THE SUPPORT I NEED AND SO MUCH MORE! I LOVE HOW EASY IT IS TO ADJUST MY POSITION!

– EMILY R.

HIGH-RESILIENCY FIBER FILL  
for maximum comfort

FABRIC  
Soft, comfortable, luxurious fabrics are durable as they are stylish. Tested to commercial standards, guaranteed to last.



UNIQUE SPRING DESIGN  
Interlocking springs ensure maximum support while delivering individualized comfort.


STEEL FRAME  
Heavy-duty metal framework is forged to last.



STRONG WOOD FRAME  
Premium furniture-grade hardwood is robust, durable and long-lasting.

MORTISE TENON JOINTS  
Interlocking puzzle-like joints

HEAVY-DUTY LIFT ACTUATOR AND SCISSOR MECHANISM



Design also includes integrated battery backup and quiet and smooth lift system.



You deserve to have a room that's entirely your own. Fill it with your favorite mementos, stylish décor and most importantly, a comfortable and luxurious lift chair. The Oasis collection by Pride® Power Lift Recliners delivers the ultimate in relaxation with an overstuffed, biscuit-back design and contoured memory foam seating. Experience true-infinite positioning and program your perfect position with the memory button hand control. Enjoy a refreshing beverage and savor the moments of absolute solitude.

Feel the luxury of beautiful Saratoga fabrics under your fingers and find the perfect match. With convenient dual pockets for storage and a standard lumbar pillow for additional support, the Oasis collection brings a whole new level of comfort to your home. It's time to immerse yourself in your oasis with Pride Power Lift Recliners.

# OASIS



## STANDARD FABRIC: SARATOGA



CASHMERE CHARCOAL EGGPLANT GODIVA NAVY



When it's been a long day and you need to rejuvenate your body, mind and spirit, consider the power of true-infinite positioning. The Infinity collection by Pride® Power Lift Recliners is designed to be fully adjustable to your comfort needs. At the touch of a button on the sleek hand remote, the backrest and footrest move independently, so you get the support you need wherever you need it.

With pocketed coils that conform to your body, feel the stress of the day fade away. The Infinity collection is the ideal addition to your living room or family room. Available in four different fabric selections, the Infinity Recliner is built to last, year after year. As your life changes, one thing remains the same. The absolute comfort and support of an Infinity Power Lift Recliner.



STANDARD FABRIC: DURASOFT

# INFINITY



TIMBER

OAT

DEEP SKY

EMBER



There is a chill in the morning air as a new day begins. Before taking on the bustle of your life, embrace the inviting warmth and coziness of the Heritage collection by Pride® Power Lift Recliners. Stretch out completely with the footrest extension and enjoy the steaming aroma from your morning coffee. Total Comfort Seating ensures outstanding relaxation and comfort, as you feel the softness of luxurious Cloud 9 fabrics enveloping your body.

When you can't put the day off any longer, innovative First Up® technology on the Heritage collection lifts you quickly and smoothly, helping you to get up and get moving. With a quiet and smooth lift system that offers durable and long-lasting performance, the Heritage collection fulfills your every need. For the recliner that perfectly complements every aspect of your life, choose the Heritage collection by Pride.



STANDARD FABRIC: CLOUD 9

# HERITAGE





## LC-105

When it comes to the essentials of everyday life, never sacrifice your comfort. Experience a simply elegant solution with the Essential LC-105 Power Lift Recliner. Fill the empty spot in your house or apartment with the chair that complements your décor perfectly. With beautiful micro-suede fabrics in soft, dark colors, the Essential Lift Chair from Pride® Power Lift Recliners is the accent piece you've been waiting for!

The Essential LC-105 has an easy-to-use hand remote, so finding a comfortable position is as simple as pressing a button! Sinuous Spring Seating gives you endless support. Rolled, pleated armrests add a touch of elegance. With dual pockets that are easily accessible, the essentials of your life are always at your fingertips. Discover the elegant simplicity of the Essential LC-105 Recliner.



STANDARD FABRIC: MICRO-SUEDE




COCOA



SANDAL



SKY



NEW CHESTNUT

OPTIONAL FABRIC:  
LEXIS URETHANE®



# ESSENTIAL



## LC-250

Although the demands of life are relentless, take solace in the comfort of a Pride® Power Lift Recliner. The Essential LC-250 Recliner gives you the boost you need to take on the day. Total Comfort Seating provides complete support for ultimate relaxation. The two-button hand remote is simple to use. Adjust your chair effortlessly to find your perfect position every time.


With luxurious Cloud 9 fabrics that are soft to the touch, you'll find the perfect match to suit your home décor. Enjoy the plentiful storage of dual pockets that are easy to reach. Experience quiet and smooth lift performance, day after day, year after year. It's time to invest in your comfort with the Essential LC-250 Recliner.



STANDARD FABRIC: CLOUD 9



BLACK CHERRY



PACIFIC



STONE



WALNUT



NEW CHESTNUT

OPTIONAL FABRIC:  
LEXIS STA-KLEEN®

# ESSENTIAL

ULTRAFABRICS® ULTRALEATHER®



Ultrafabrics® Ultraleather® on Pride® Power Lift Recliners offer a richly detailed grain pattern expressive of premium grade leather. Ultrafabrics polyurethane products are manufactured utilizing Tukumi™ technology, a 4-layer laboratory process to create a luxurious, durable and comfortable polyurethane fabric.

FEATURES AND BENEFITS

- Rich feel and natural grain pattern expressive of soft luxury-grade leather
- Clinically proven more resistant to temperature changes than vinyl and genuine leather for comfort
- Antimicrobial surface cleans effortlessly with soap and water

CRYPTON® SUPER FABRIC



Crypton® Super Fabrics from Pride® Power Lift Recliners are not only luxurious, durable and contemporary, but also resist spills, stains, odors and bacteria. Nothing gets past Crypton's unique patented integrated fiber process, protecting your upholstery inside and out.

FEATURES AND BENEFITS

- 100% water resistant – spills just wipe away
- Only disinfectable fabrics in the U.S. approved by the EPA
- Permanent technology that will not come off the fabrics
- 5-year warranty on moisture barrier

LEXIS STA-KLEEN® PERFORMANCE FABRIC



Lexis Sta-Kleen® Performance Fabric is a revolution in faux leather upholstery — the first stain resistant urethane fabric. Sta-Kleen offers permanent protection, created with a proprietary bonding process inherent in the manufacturing process of the product. It is not a topical application; it will last the lifetime of the upholstery.

FEATURES AND BENEFITS

- Look and feel of real leather, but with care-free durability
- Exceptional strength, abrasion resistance, and UV stability
- “Stain shield” allows spills to wipe off with a cloth, and it's anti-bacterial

QUICK REFERENCE FABRIC GUIDE

MODEL	SARATOGA	DURASOFT	CLOUD 9	MICRO-SUEDE	LEXIS STA-KLEEN	CRYPTON ARIA	ULTRALEATHER	LEXIS URETHANE
LC-580iM	Standard	N/A	N/A	N/A	Optional	Optional	Optional	N/A
LC-580iL	Standard	N/A	N/A	N/A	Optional	Optional	Optional	N/A
LC-525iPW	N/A	Standard	N/A	N/A	Optional	Optional	Optional	N/A
LC-525iS	N/A	Standard	N/A	N/A	Optional	Optional	Optional	N/A
LC-525iM	N/A	Standard	N/A	N/A	Optional	Optional	Optional	N/A
LC-525iL	N/A	Standard	N/A	N/A	Optional	Optional	Optional	N/A
LC-358PW	N/A	N/A	Standard	N/A	Optional	Optional	Optional	N/A
LC-358S	N/A	N/A	Standard	N/A	Optional	Optional	Optional	N/A
LC-358M	N/A	N/A	Standard	N/A	Optional	Optional	Optional	N/A
LC-358L	N/A	N/A	Standard	N/A	Optional	Optional	Optional	N/A
LC-358XL	N/A	N/A	Standard	N/A	Optional	Optional	Optional	N/A
LC-358XXL	N/A	N/A	Standard (only Stone & Walnut)	N/A	Optional (only Mushroom & Chestnut)	N/A	N/A	N/A
LC-105	N/A	N/A	N/A	Standard	N/A	N/A	N/A	Optional
LC-250	N/A	N/A	Standard	N/A	Optional (only Chestnut)	N/A	N/A	N/A

ACCESSORIES



**NON-SKID LEG LEVELERS**  
Non-skid leg levelers keep your Power Lift Recliner from sliding on a vast array of flooring.



**HEAD & LUMBAR PILLOW (SOLD SEPARATELY)<sup>4</sup>**  
Adds additional support and comfort (in the same fabric the chair was ordered in) in the head or lumbar area.



**DELUXE HEAT AND MASSAGE/HEAT AND MASSAGE**  
Deluxe heat and massage features four massagers and two heat pads. Regular heat and massage offers two massage units and one heat pad.



**FIRST UP® TECHNOLOGY<sup>1</sup>**  
Patent-pending technology offers up to a 50% increase in lift speed vs conventional lift chairs.



**CLEANING KIT**  
This kit can be used with any Power Lift Recliner fabric (excluding vinyl/leather) in order to remove almost any stain or odor.



**WIRELESS REMOTE<sup>2</sup>**  
Performs lifted, seated, tv and full recline chair positions.

COLLECTION	OASIS	OASIS	INFINITY	INFINITY	INFINITY	INFINITY
MODEL	LC-580iM	LC-580iL	LC-525iPW	LC-525iS	LC-525iM	LC-525iL
POSITION	True-Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'4" - 6'	True-Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'9" - 6'2"	True-Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'2" & Below	True-Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'3" & Below	True-Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'4" - 6'	True-Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'9" - 6'2"
SIZE	Medium	Large	Petite Wide	Small	Medium	Large
WEIGHT CAPACITY	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.
SEAT-TO-FLOOR	20"	21"	19"	19"	20"	21"
SEAT DEPTH	21"	21"	19"	19"	20"	22"
SEAT WIDTH	20"	21"	24"	19"	20"	22"
TOP OF BACK TO SEAT	26"	29"	28"	26"	28"	30"
FULLY PADDED CHAISE	Yes	Yes	Yes	Yes	Yes	Yes
BACK STYLE	Blown Fiber Biscuit Back	Blown Fiber Biscuit Back	Blown Fiber Split Back	Blown Fiber Split Back	Blown Fiber Split Back	Blown Fiber Split Back
ADDITIONAL OPTIONS	Deluxe Heat & Massage <sup>3</sup>	Deluxe Heat & Massage <sup>3</sup>	Deluxe Heat & Massage <sup>3</sup>	Deluxe Heat & Massage <sup>3</sup>	Deluxe Heat & Massage <sup>3</sup>	Deluxe Heat & Massage <sup>3</sup>
FOOTREST EXTENSION	N/A	N/A	N/A	N/A	N/A	N/A
SEATING TYPE	Pocketed Coil Springs w/ Memory Foam Topper	Pocketed Coil Springs w/ Memory Foam Topper	Pocketed Coil Springs	Pocketed Coil Springs	Pocketed Coil Springs	Pocketed Coil Springs
DISTANCE FROM WALL	20"	21"	19"	20"	20"	21"
OVERALL WIDTH	39"	39"	38"	34"	34.5"	37"
CHAIR WEIGHT	130 lbs.	132 lbs.	130 lbs.	118 lbs.	130 lbs.	130 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated


\*Not all units achieve all positions. Please refer to specific product.

<sup>1</sup> Available on Heritage Collection except LC-358XL & XXL | <sup>2</sup> Available on LC-358 S, PW, M & L | <sup>3</sup> USB hand control not available with heat & massage options. | <sup>4</sup> Not available on all models.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

HERITAGE	HERITAGE	HERITAGE	HERITAGE	HERITAGE	HERITAGE	ESSENTIAL	ESSENTIAL
LC-358PW	LC-358S	LC-358M	LC-358L	LC-358XL	LC-358XXL	LC-105	LC-250
3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'4" - 6'	3-Position Chaise Lounger 5'9" - 6'2"	3-Position Chaise Lounger 5'6" - 6'1"	2-Position Chaise Lounger 5'6" - 6'1"	3-Position Chaise Lounger 5'2" - 5'10"	3-Position Chaise Lounger 5'4" - 6'
Petite Wide	Small	Medium	Large	Heavy Duty	Very Heavy Duty	Medium	Medium
375 lbs.	375 lbs.	375 lbs.	375 lbs.	500 lbs.	600 lbs.	325 lbs.	375 lbs.
19"	19"	20"	21"	20"	20"	18"	20"
19"	19"	20"	22"	22"	22"	19"	20"
24"	18"	20"	22"	26"	30"	20"	22"
23.5"	24"	24.5"	25.5"	29"	29"	28"	26"
Yes	Yes	Yes	Yes	Yes (removable)	Yes (removable)	Yes	Yes
Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Sewn Pillow Back	Blown Fiber Split Back
Heat & Massage <sup>3</sup> , First Up <sup>®</sup> , Wireless Remote	Heat & Massage <sup>3</sup> , First Up <sup>®</sup> , Wireless Remote	Heat & Massage <sup>3</sup> , First Up <sup>®</sup> , Wireless Remote	Heat & Massage <sup>3</sup> , First Up <sup>®</sup> , Wireless Remote	Heat & Massage <sup>3</sup>	No	No	No
N/A	Standard	Standard	Standard	N/A	N/A	N/A	N/A
Total Comfort	Total Comfort	Total Comfort	Total Comfort	Sinuous Spring	Sinuous Spring	Sinuous Spring	Total Comfort
18"	18"	17"	19"	19"	14"	17"	20"
40"	35"	35"	37"	39.5"	43"	33"	34"
130 lbs.	118 lbs.	117 lbs.	120 lbs.	159 lbs.	168 lbs.	96 lbs.	118 lbs.
7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	1-Yr Limited	2-Yr. Limited

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.



PRIDE<sup>®</sup>  
POWER  
RECLINER  
COLLECTION

(US) 800-800-8586 | (CAN) 888-570-1113 | [pridemobility.com](http://pridemobility.com) | [prideliftchairs.com](http://prideliftchairs.com) |    

