

# SOFTCAIRE SIMPLE DIGITAL W/ EC MATTRESS

Weight Capacity: 350 lbs



## Alternating Pressure with Cell on Foam Design Support

Our Best-in-Class Offering, the "Softcaire with Simple Digital and EC Mattress" includes 17, 5" alternating air cells made of durable PVC, with Cell-on-Foam Design (3" Foam Base) for Effective Pressure Redistribution, Power Outage Protection and Patient Transport

## Simplified Digital Pump

8 lpm Air Output, Auto-Sensing Tech, and **Static Function with Auto-Return** so that you can "Set it and Forget it!"

## Cover

2-Way Stretch Top Cover that is Low Shear, Water-Resistant, Vapor-Permeable, and Anti-Microbial

## Quick Connector Cap

Ensures mattress remains inflated during power outage or patient transport

## CPR Function

Allows for one-handed, rapid deflation

## Pillow Function

Head zone comprised of three fully inflated static air cells to provide patient comfort

## Spring-Loaded, Anti-Kink Tubes

Enhances experience and reduces service calls

### Pump Specifications

Pump:	Digital Control
Air Output:	8 liter/min
Pressure Range:	20-55 mmHg
Cycle Time:	10 min
Power Supply:	120V/60 Hz

### Mattress Specifications

Inflated Dimension:	80" x 36" x 5" + 3" Foam Base
Alternating:	A/B Alternating
Structure:	17 x 5" Air Cells + 3" Foam Base
Cover Material:	Nylon/PU Quilted with Zipper
Cell & Base Material:	Nylon/PVC
Weight Capacity:	350 lbs
HCPCS	E0277

**SET IT AND FORGET IT!**  
With Static Function with Auto-Return



# SOFTCAIRE SIMPLE DIGITAL W/ EC MATTRESS

Weight Capacity: 350 lbs

## 18 MONTH LIMITED WARRANTY

Your Rhythm product is warranted to be free of defects in materials and workmanship for 18 Months to the original consumer purchaser. This warranty does not cover device failure due to misuse, negligence, improper storage or handling, improper operation and unauthorized modifications. This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, which are subject to normal wear and need periodic replacement.

If the product is to be returned to Rhythm Healthcare, prior authorization will be required. This product is identified by a serial number at the time of manufacture. Removing or defacing the serial number label may void this warranty. If you have a question about your Rhythm Healthcare product or this warranty, please email us at [contactus@rhythmhc.com](mailto:contactus@rhythmhc.com)