

SOFTCAIRE DX BARIATRIC (DELUXE SYSTEM)

Weight Capacity: 660 lbs



ITEM #	DESCRIPTION
H8841142	SoftCaire DX Bariatric (Deluxe System)

Alternating Pressure with Cell on Cell Design Support

- Our Premium Offering, the Dual-Therapy “Delux Bariatric” includes 20, 10” alternating air cells made of durable TPU, with Laser Drilled Low Air Loss Therapy Embedded Cell-on-cell Design for Effective Pressure Redistribution, Power Outage Protection and Patient Transport

Premium Digital Pump

- 12 lpm Air Output
- Auto-Sensing Tech Adjusts Air Pressure to Patient’s Weight
- Adjustable Cycle Time, Seat-Inflate, Static and Max-Firm Modes

Cover

- 2-Way Stretch, Top Cover is Low Shear, Water-Resistant, Vapor-Permeable, and Anti-Microbial

Quick Connector with Cap

- Ensures Mattress Remains Inflated during power outage or patient transport

CPR Function

- Allows for One-Handed, Rapid Deflation in case of medical emergencies or to facilitate rapid disassembly

Pillow Function

- Head Zone Comprised of Three Fully Inflated Static Air Cells to provide patient comfort

Heel Function

- Provides Single Cell Deflation to Float Vulnerable Heels

Spring-Loaded, Anti-Kink Tubes

- Enhances experience and reduces service calls

SOFTCAIRE DX BARIATRIC (DELUXE SYSTEM)

Weight Capacity: 660 lbs



SPECIFICATIONS

Pump

Weight Setting Control:	Digital
Pump Dimension:	13.4"L x 5"W x 9.6"H
Weight:	7.75 lbs
Air Output:	12 liter/min
Case Material:	Flame Retardant Abs
Pressure Range:	20-65 mmHg
Cycle Time:	10/15/20/25 min
Power Supply:	120V/60 Hz
Power Cord Length:	14' 9"

Mattress

Inflated Dimension:	80"L x 42"W x 10"H
Alternating:	A/B Alternating
Structure:	20 x 10" Air Bladders
Cover Material:	Nylon/PU Quilted with Zipper
Cell Material:	Nylon/TPU
Base Material:	Durable Nylon/PVC 840D
Maximum Weight Capacity:	660 lbs
Flame Retardant Standards:	CA117

18 MONTH LIMITED WARRANTY

Your Rhythm product is warranted to be free of defects in materials and workmanship for eighteen (18) months to the original consumer purchaser. This warranty does not cover device failure due to misuse, negligence, improper storage or handling, improper operation and unauthorized modifications. This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, which are subject to normal wear and need periodic replacement.

If the product is to be returned to Rhythm Healthcare, prior authorization will be required. This product is identified by a serial number at the time of manufacture. Removing or defacing the serial number label may void this warranty. If you have a question about your Rhythm Healthcare product or this warranty, please email us at contactus@rhythmhc.com