

Scout Compact 3-Wheel Scooter

HCPCS: K0800-1,K0800



FIGURE A




FIGURE B



FIGURE C




FIGURE D




FIGURE E



FIGURE F

The Scout transportable scooter combines easy and convenient mobility with superior value. The Scout is compact, lightweight and easy to disassemble-perfect for taking wherever you want to go, either around the neighborhood or far from home. It can be disassembled quickly and easily into five easy-to-handle pieces for loading into a car, truck, or van and taking shopping or on vacation. The three-wheel version of the Scout allows for nimbly negotiating indoor obstacles. The Scout also offers great features like a quick connect battery pack and a delta-style tiller to make operating the scooter easier for those with limited strength or dexterity. This model comes with color-through panels to hide scuffs and scratches and black wheel rims.

- Quick and easy disassembly
- Interchangeable color panels hide scuffs and scratches
- Armrests are padded and adjustable
- Quick connect batteries
- Ergonomic throttle control
- Height-adjustable swivel seat with fold-down backrest
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust tiller
- Anti-tip wheels
- Delta tiller
- Easy freewheel operation

PRODUCTS	Description
SFSCOUT3	Scout 3-Wheel Scooter, Red/Blue, 1/ea
SPECIFICATIONS	Description
Warranty	24 Month Limited Warranty on Electronics
Product Weight	85 lb
Product Dimensions	42.5" (L) x 20.5" (W) x 34" (H)
Seat Dimensions	16" (D) x 16.5" (W)
Weight Capacity	300 lb
Top Speed	4.25 mph
Batteries (pair)	12V x 12AH
Battery Charger	2A Offboard
Climbing Angle	6°
Color	Red & Blue
Controller	Dynamic 45A
Ground Clearance	2.5"

SPECIFICATIONS	Description
HCPCS	K0800
Maximum Range	9 miles
Motor	24V x 270W x 4700 RPM
Seat Type	Stadium
Turning Radius	45.5"
Warranty on Batteries	12 Month Warranty on Batteries
Warranty on Frame	Limited Lifetime Warranty on Frame
Wheel Size	Front: 8" x 2" Rear: 8" x 2"
ACCESSORIES	Description
SF8030	Crutch Holder Set, SType, 1/ea, EA
PARTS	Description
SC31053	Transaxle W/ Motor & Brake, Scout, 1/ea, EA
SC31052	Controller, 45A, Scout 3-Wheel, 1/ea, EA
SC31042-D	Blue Decorations, Scout 3 Deluxe, 1/ea, EA
SC31040	Front Frame, Scout 3, 1/ea, EA
SC31041-D	Red Decorations, Scout 3 Deluxe, 1/ea, EA
SC31042	Blue Decorations, Scout 3, 1/ea, EA
SC31039	Front Shroud, Scout 3, 1/ea, EA
SC31041	Red Decorations, Scout 3, 1/ea, EA
SC31044-D	Front Wheel, Scout 3 Deluxe, 1/ea, EA
SC31046	Front Axle, Scout 3, 1/ea, EA
SC31038	Front Fork, Scout 3, 1/ea, EA
SC31043	Front Wheel, Black, Scout 3, 1/ea, EA
SC31047	Battery Box Hold-Down, Scout 3 & 4, 1/e, EA
SC31049	Rear Axle Bolt Assmbly, Scout 3 & 4, 1e, EA
SC31033	Conn.blckSet,Rear&Frntw/harn,Scout,1/ea, EA
SC31026	WireHarness,Scout,CurtisElectronics,1/ea, EA
SC31012	Rear Frame, Scout, 1/ea, EA
SC31025	Hub Cap, Scout, 1/ea, EA
SC31011	Rear Shroud, Scout, 1/ea, EA
SC31037	Electromagnetic Brake, Scout, 1/ea, EA