



For easy ordering and information, please
contact Customer Care Toll Free at
1(800) 743-0772 Fax: 1(877) 797-7081

Streamer Sport

Name: _____
Phone Number: _____
Acct. Number: _____
PO Number: _____
Address: _____
City: _____ State: _____ Zip Code: _____

Ship to: _____
Phone Number: _____
Address: _____
City: _____ State: _____ Zip Code: _____

Base Scooter	Description	MSRP	
☐ 888WA-RED	Streamer WA (prejudgment)	K0823	\$5,295.00
	Streamer WA w/ Pan Seat Upgrade	K0822	
BATTERY			
BAT-U1	Battery U1	x2	Incl.
SEATS (W"xD"xH")			
SV191821B	19X18 Captain Seat BLK w/o Recline		Incl.
ARM REST ASSY			
104300-00106-A3	(L) Flipback Height Adj. Armrest		Incl.
104400-00106-A3	(R) Flipback Height Adj. Armrest		Incl.
DRIVE TIRES			
TSR26085 Rear	Grey 10" Flat Free Drive Tires	x2	Incl.
CASTORS			
1351-00601-01(K)	6" Castor Assy w/fork Left		Incl.
1351-00602-01(K)	6" Castor Assy w/fork Right		Incl.
FOOT REST			
2329-00100-11	Flip Up Foot Plate		Incl.
CONTROLLER			
8301-00200-06	Shark Joystick Controller		Incl.
1361-00204-04	Shark Main Control Box		Incl.
CHARGER			
109201-88124	5 Amp Charger Off-Board		Incl.
ACCESSORIES			
300909-03	Safety Belt		Incl.
1329-00700-02	Joystick Knob		Incl.
Total			_____

SPECIFICATIONS

OVERALL DIMENSIONS

Base Length 36"
Base Width 24"
Height 45"

Ground to Seat Top

Ground Clearance 2"

WEIGHTS

User Weight Capacity 300 lbs

Total Unit Weight 203 lbs

Battery Weight Module 23 lbs

Tires

Number, Size of Front Tire 2, 6"

Number, Size of Rear Tire 2, 10"

Batteries

BAT - U1 Module Capacity 12v32Ahx2

Charger

Charger-On Board 5 amp

Motor Power Output

Horsepower60 hp x2

Driving Specs

Turning Radius 17.7"

Range 10 miles

Maximum Speed 3.75 mph

Maximum Climbing Angle 10% / 6°



Replacement Parts & Accessories

Part Number	Description		MSRP
BATTERY			
BAT-U1	Battery U1	x2	\$140.00ea
SEATS (W"xD"xH")			
SV191821B	19X18 Captain Seat BLK w/o Recline		\$280.00
☐ SV191820RG	19x18 Recline High-Back Grey Vinyl		\$400.00
☐ SV201820RB	20x18x20 Recline Blk Pillow Top		\$370.00
☐ SV2121211/2RG	21x21x211/2 Recline Capt Seat-Grey Vinyl		\$400.00
SV2121211/2RB	21x21x211/2 Recline Capt Seat-Blk Vinyl		\$400.00
☐ SP181821RG	18x18x21 Pan Seat;2-tone Grey		\$400.00
☐ SP202023RG	20x20x23 Pan Seat;2-tone Grey		\$400.00
☐ SP161618RG	16x16x18 Pan Seat;2-tone Grey Recline		\$400.00
ARM REST ASSY			
☐ 104300-00106-A3	(L) Flipback Height Adj. Armrest		\$140.00
☐ 104400-00106-A3	(R) Flipback Height Adj. Armrest		\$140.00
☐ 104310-00102	Arm pad w/ height-adjustment tube	x2	?
DRIVE TIRES			
☐ TSR26085 Rear	Grey 10" Flat Free Drive Tires	x2	\$50.00ea
☐ TAR26085 Rear	Grey 10" Pnuematic Tires	x2	\$32.00ea
☐ 105102-88100	10" Tire Tubes	x2	\$20.00ea
CASTORS			
☐ 1351-00601-01(K)	6" Castor Assy w/fork Left	x2	?
☐ 1351-00602-01(K)	6" Castor Assy w/fork Right	x2	?
☐ 205201-00600-A2	6" Castor Tire Only	x2	\$90.00ea
FOOT REST			
☐ 2329-00100-11	Flip Up Foot Plate		\$150.00
☐ ELR-888WA-Red	ELR Kit complete Streamer-888WA		\$180.00
CONTROLLER			
☐ 8301-00200-06	Shark Joystick Controller		\$360.00
☐ 1361-00204-04	Shark Main Control Box		\$776.00
CHARGER			
☐ 109201-88124	5 Amp Charger Off-Board		\$230.00
ACCESSORIES			
☐ BKS100L(K)	Left Stump Support		\$340.00
BKS100R(K)	Right Stump Support		\$340.00
☐ 1343-00709-01-A3	Right Swing Away Mount Bracket for VSI		\$80.00
☐ 1343-00709-00-A3	Left Swing Away Mount Bracket for VSI		\$80.00
☐ O210(K)	Oxygen Tank Holder (Metal)		\$150.00
WH-0185(K)	Walker Holder		\$150.00
☐ OXY-03(K)	Mesh Oxygen Tank Holder		\$60.00
☐ 300909-03	Safety Belt		\$34.00
☐ 1329-00700-02	Joystick Knob		\$12.00

Total