

VivaLIFT![®]
POWER RECLINERS
By Pride® Mobility Products Corp.

THE MOST
COMFORTABLE
RECLINER
IN THE WORLD!

THE 4-MOTOR RECLINER

WITH
COUNTLESS
POSITION
COMBINATIONS



Shown in Saville
Brown Fabric

METRO

Contemporary comfort has never been more stylish! The Metro Collection by **VivaLift!**[®] Power Recliners offers maximum relaxation with infinite, lay-flat positioning. Soft, dark colors with Saville fabric offers a unique and cozy appeal, making the Metro recliner a perfect accent piece for the bright areas in your home. A USB port in the remote offers convenient access for charging smart devices. It's time to experience the richness and luxury of the VivaLift![®] brand.



Pride Mobility[®] (US) 800-800-8586 | (CAN) 888-570-1113
pridemobility.com | prideliftchairs.com |    

METRO

COLLECTION	METRO
MODEL	PLR-925M
POSITION	Infinite Lay Flat Relaxed Reclined Upward Lift
SIZE	Medium
WEIGHT CAPACITY	375 lbs.
SEAT-TO-FLOOR	19"
SEAT DEPTH	21"
SEAT WIDTH	21"
TOP OF BACK TO SEAT	26.5"
FULLY PADDED CHAISE	Yes
FOOTREST EXTENSION	Standard
SEATING TYPE	Total Comfort
CHAIR HEIGHT	41"
DISTANCE FROM WALL	14"
OVERALL WIDTH	33"
CHAIR WEIGHT	142 lbs.
WARRANTY	7-Yr. Pro-rated

FUNCTIONS AND FEATURES:

-  POWER HEADREST
-  POWER LUMBAR
-  USB REMOTE
-  LITHIUM BATTERY BACKUP
-  FOOTREST EXTENSION
-  POWER BACKREST

POSITION GUIDE:

-  RELAXED
-  RECLINED
-  UPWARD LIFT
-  INFINITE LAY FLAT

STANDARD FABRIC: SAVILLE



GREY

BROWN



Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

 = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

Pride Mobility® (US) 800-800-8586 | (CAN) 888-570-1113

pridemobility.com | prideliftchairs.com |    