



S542

Portable Folding Electric Scooter

S542 Portable Folding Electric Scooter

Features

- Strong light-weight alloy aluminum frame for easy folding and carrying.
- Small foldable dimension for easy storage.
- Easy attach 24V Li-Ion Battery, Quick charge with 2A off-board charger.
- Height and angle adjustable steering tiller to fit individual user demand.
- Ideal for both indoor and outdoor use.
- Patented Twin-front-wheels w/spring suspension system offers excellent stability during sharp turning.
- Optional Armrest.


Dimension (cm / in)

| | | |
|-----------------------|-------------|-------------------------|
| | | |
| seat width | seat depth | seat height from ground |
| 40.5/16 | 35/13.8 | 59.5/23.4 |
| | | |
| seat height from deck | back height | wheel base |
| 39/15.4 | 31.5/12.4 | 78.5/31 |

| | | |
|------------------------------------|---------|---------|
| | | |
| length | width | height |
| Unfold | 98/38.6 | 43.5/17 |
| Folded | 52/20.5 | 43.5/17 |
| tiller: 102.5/40.3 ; seat: 87/34.2 | | |
| 64/25.2 | | |

Weight (kg / lb)

| base w/seat (w/o batt. and armrest) | battery/ea | ttl w/batt. (w/o armrest) | capacity |
|--|------------|------------------------------|----------|
| 24/53 | 1.7/3.7 | 25.7/56.7 | 113/250 |



Technical Specs.

| | |
|------------------|-----------------------------------|
| max. speed up to | 6 kph / 4mph |
| range up to | 18km / 11.25 mile |
| turn radius | 96cm / 37.8" |
| ground clearance | 9.5cm / 3.7" |
| motor | DC24V, 80W (Rating) / 340W (Max.) |
| controller | PG S-Drive 45 |
| battery | Li-Ion Battery 24V / 11.5AH |
| charger | 2A off-board |
| gradient | 6° |
| front wheel | 6" (153*40) foam filled tire |
| rear wheel | 7" (178*45) foam filled tire |
| brake | Electromagnetic Brake |



No. 18, Jingke Road, T.P.M.T. Park, Nantun District,
Taichung City 40852, Taiwan, R.O.C.
TEL: +886-4-23594991 / FAX: +886-4-23594992
www.merits.com.tw / sales@merits.com.tw

