

Junior (P320) Group 2
 Super Light Mini RWD Power Wheelchair

Order Form (Effective April 19, 2017)

For ease of ordering contact Inside Sales toll free at 1-800-963-7487, Fax: 239-574-2661



Order Taken by: _____

Company Name: _____ Account # _____

Phone: _____ Date: _____ P.O # _____

Ship to Name/Address: _____

Special Note:

DME suppliers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.

CHAIR TYPE	PRICE	HCPCS	PRODUCT SPECIFICATIONS
Junior Group 2			
<input type="checkbox"/> P320-2S C Rear Wheel Drive Base (18"W x 16"D)	\$6275.00	K0823	Dimensions Length 33" Width 24" Height 35" Ground Clearance 2" Seat Dimensions Width 18" Depth 16" Back Height 17" Weight Capacity 300-lbs Maximum Speed 4.5mph Range Up To 12.5 miles Tires Front/Rear Casters 7 x 1 3/4" Drive Wheel 9" Turning Radius 23" Incline Capability 8° Product Weight Total Weight without Battery 105.5-lbs Total Weight with Battery 135-lbs Weight of Heaviest Piece 65.5-lbs Battery Type 12V / 22A gel x 2 Weight 29.5-lbs x 2 Warranty Frame 5 year / Limited Lifetime Electronics 1 year Motor/Brakes 1 year Batteries 6 months
<input type="checkbox"/> P320-2S S Rear Wheel Drive Base (18"W x 16"D).	\$6475.00	K0822	
FRONT RIGGINGS			
<input type="checkbox"/> Footplate*	No Charge		
<input type="checkbox"/> Elevating Leg Rests.....	\$360.00pr	E0990	
JOYSTICK MOUNTING HARDWARE			
<input type="checkbox"/> Swing Away Joystick Mount	\$240.00	E1028	
<input type="checkbox"/> Right Mount <input type="checkbox"/> Left Mount.....			
BATTERY OPTIONS			
<input type="checkbox"/> 12V / 22A Battery (On Chair).....	\$85.00ea	K0733	
ARM RESTS			
<input type="checkbox"/> Flip Up, Adjustable-Height	\$150.00ea	E0973	
SEATING OPTIONS (for P320-2S S)			
<input type="checkbox"/> 18"Wx16"D Promed Gen.Use Seat Cushion 71002	\$127.92	E2601	
<input type="checkbox"/> 18"W x 16"D Promed Gel Seat Cushion #74002	\$199.92	E2603	
OPTIONS & ACCESSORIES			
<input type="checkbox"/> Oxygen Tank Holder.....	\$138.00	E2208	
<input type="checkbox"/> Cane & Crutch Holder.....	\$98.00	E2207	
<input type="checkbox"/> Walker Holder.....	\$118.00	K0108	
<input type="checkbox"/> Foldable Rear Basket (Med.).....	\$72.00		
<input type="checkbox"/> Foldable Rear Basket (Lrg.).....	\$108.00		
<input type="checkbox"/> Canopy.....	\$298.00		
<input type="checkbox"/> Cover.....	\$180.00		
<input type="checkbox"/> Cup Holder**.....	\$24.00		
<input type="checkbox"/> Light**.....	\$64.00		
<input type="checkbox"/> Mounting Bracket.....	\$15.00		
<input type="checkbox"/> Side Bag.....	\$18.00		

*Standard with Unit
 **Mounting Bracket Required

SHROUD COLORS*

Red..... Blue..... No Charge



Merits Health Products, Inc.,
 730 NE 19th Place,
 Cape Coral, FL 33909
www.meritshealth.com

Quantity for above selection

Features:

Light weight and easy to disassemble for handling, storage, and transportation
 Programmable controller with a wide range of selections to suit the needs and demands of different users
 Solid Caster Wheels & Drive Wheels
 3A off-board charger
 Lap Belt (Standard with chair)