

SHOPRIDER®
A TRUSTED NAME IN MOBILITY



Scooter Catalog

Echo Folding

The Echo Folding is designed for the user on the go. The scooter folds for transportation and unfolds for use in a matter of seconds. The unit weighs in at just 59 lb. and can store easily in small spaces. It comes standard with solid tires, electric brake system and a curved tiller for added comfort. Available in red, gray and black, the Echo Folding will get you to your destination in style.

● Red
 ● Gray
 ● Black

Model Number: FS777



The Echo Folding easily folds down for storage or transportation.

Technical Specifications

Overall Dimensions (L x W x H)	in	39.5" x 17.5" x 35"
Number, Size of Tires	Front, Rear	2, 6" (2, 8")
Suggested User Weight ¹	lb	Max. 200
Weight of Heaviest Part ²	lb	59
Weight ³	lb	59
Each Battery Module Weight	lb	7
Battery Size	Type x PCS	12V10AH x 2
Charger	Amp	Off board (2A)
Output Power of Motor	HP	0.45
Maximum Speed	MPH	3.75
Maximum Climbing Angle/ Safe Climbing Angle ⁵	%(°)	10(6)
Range (Per Charge) ⁶	Mile	8.7-10
Turning Radius	in	36"
Ground Clearance	in	1.5"
Handle Bar	Type	T-Bar
Model #		FS777
Standard Seat	W x D x BH	SF151413B

Echo

The Echo scooter has a light weight frame designed for the user on the go. The scooter disassembles into three pieces for portability and storage; with the heaviest part weighing just 35 lb. It has a height adjustable seat, allowing for improved comfort during driving operation. The Echo 3 also comes standard with a quick connect battery pack, solid tires and a basket to store personal items.

- + 250lb weight capacity
- + Quick connect battery pack
- + Removable battery for charging
- + Adjustable seat height



● Red
 ● Blue

Model Number: SL73

Technical Specifications

Overall Dimensions (L x W x H)	in	37" x 21" x 36"
Number, Size of Tires	Front, Rear	1, 8" (2, 8")
Suggested User Weight ¹	lb	Max. 250
Weight of Heaviest Part ²	lb	35
Weight ³	lb	70
Each Battery Module Weight	lb	15
Battery Size	Type x PCS	12V10AH x 2
Charger	Amp	Off board (2A)
Output Power of Motor	HP	0.45
Maximum Speed	MPH	3.8
Maximum Climbing Angle/ Safe Climbing Angle ⁴	%(°)	10(6)
Range (Per Charge) ⁶	Mile	8
Turning Radius	in	32"
Ground Clearance	in	1.5"
Handle Bar	Type	T-Bar
Model #		SL73
Standard Seat	W x D x BH	SF171413B

Note: ¹ On level Road. ² Without basket and seating. ³ With basket, seat & batteries. ⁴ Base on user weight of 200LB. ⁵ Based on user weight of 165LB. ⁶ After the batteries and mechanical moving parts fully break in.



Hero

The Hero scooter has a light weight frame designed for the user on the go. The scooter disassembles into three main parts for portability and storage; with the heaviest part weighing just 40 lb. It has a height adjustable seat, allowing for improved comfort during driving operation. The Hero also comes standard with a quick connect battery pack, solid tires and a basket to store personal items.

- + 250lb Weight Capacity
- + Keyed Ignition Security
- + Quick Connect Battery Pack
- + Curved Tiller
- + Swivel Seating



● Red ● Blue
 Model Number: SL73N



Scootie

The Scootie is light-weight, easy to transport and has articulating front wheels. The Scootie weighs in at just 82 lb. with the batteries installed. The patented articulating front wheels allow four-wheel stability with three-wheel maneuverability. The Scootie is length adjustable and can fold down for easy transport and storage.

- + 250lb weight capacity
- + Length Adjustable
- + Removable Battery Pack
- + Keyed Ignition Security
- + Adjustable Flip-up Arms



● Red ● Blue
 Model Number: TE-787NA
 Medicare Billing Code: K0800

HCPSC codes provided should not be considered as legal advice, and does not guarantee reimbursement.

Technical Specifications

Overall Dimensions (L x W x H)	in	39" x 22" x 35"
Number, Size of Tires	Front, Rear	1, 8" (2, 8")
Suggested User Weight ¹	lb	Max. 250
Weight of Heaviest Part ²	lb	40
Weight ³	lb	73
Each Battery Module Weight	lb	15
Battery Size	Type x PCS	12V10AH x 2
Charger	Amp	Off board (2A)
Output Power of Motor	HP	0.45
Maximum Speed	MPH	3.8
Maximum Climbing Angle/ Safe Climbing Angle ⁴	%(°)	10(6)
Range (Per Charge) ⁵	Mile	8
Turning Radius	in	33"
Ground Clearance	in	1.5"
Handle Bar	Type	T-Bar
Model #		SL73N
Standard Seat	W x D x BH	SF171413B

Note: ¹ On level Road. ² Without basket and seating. ³ With basket, seat & batteries. ⁴ Base on user weight of 200LB. ⁵ Based on user weight of 165LB. ⁶ After the batteries and mechanical moving parts fully break in.

Technical Specifications

Overall Dimensions (L x W x H)	in	37-40" x 17" x 34"
Number, Size of Tires	Front, Rear	2, 7" (2, 8")
Suggested User Weight ¹	lb	Max. 250
Weight of Heaviest Part ²	lb	25
Weight ³	lb	82
Each Battery Module Weight	lb	24
Battery Size	Type x PCS	12V14AH x 2
Charger	Amp	Off board (2A)
Output Power of Motor	HP	0.75
Maximum Speed	MPH	4
Maximum Climbing Angle/ Safe Climbing Angle ⁴	%(°)	10(6)
Range (Per Charge) ⁵	Mile	10
Turning Radius	in	32"
Ground Clearance	in	1.5"
Handle Bar	Type	T-Bar
Model #		TE-787NA
Standard Seat	W x D x BH	SF151413BG



Dasher 3

The Dasher 3 scooter is designed for the user on the go. It has a 250 lb. weight capacity and has a quick-disconnect system for effortless disassembly. The unit disassembles into three parts, the heaviest weighing just 26 lb. The Dasher 3 maneuvers with ease and has an adjustable tiller for added comfort during operation. The Dasher 3 also comes standard with solid tires for a smoother ride and swivel seating system.

- + 250lb Weight Capacity
- + Swivel Seating
- + Dismantles into Three Parts
- + Curved Tiller



● Burgundy

Model Number: GK83
Medicare Billing Code: K0800

HCPCS codes provided should not be considered as legal advice, and does not guarantee reimbursement.



Escape

The Escape scooter is designed for the user on the go. It has a 250 lb. weight capacity and has a quick-disconnect system for effortless disassembly. The unit disassembles into three parts, the heaviest weighing just 33 lb. The Escape maneuvers with ease and has an adjustable tiller for added comfort during operation. The Escape also comes standard with solid tires for a smoother ride and swivel seating system. The Escape also features a split battery pack design that make it convenient for travel.

- + 250lb Weight Capacity
- + Heaviest Part Weighs 33lb
- + Dismantles into Three Parts
- + Swivel Seating
- + Curved Tiller



● Red ● Blue

Model Number: 7A

Technical Specifications

Overall Dimensions (L x W x H)	in	41.5" x 21.5" x 35"
Number, Size of Tires	Front, Rear	1, 8" (2, 8")
Suggested User Weight ¹	lb	Max. 250
Weight of Heaviest Part ²	lb	25
Weight ³	lb	87
Each Battery Module Weight	lb	22
Battery Size	Type x PCS	12V14AH x 2
Charger	Amp	Off board (2A)
Output Power of Motor	HP	0.5
Maximum Speed	MPH	3.8
Maximum Climbing Angle/ Safe Climbing Angle ⁴	%(°)	10(6)
Range (Per Charge) ⁵	Mile	7.5
Turning Radius	in	33"
Ground Clearance	in	1.5"
Handle Bar	Type	T-Bar
Model #		GK83
Standard Seat	W x D x BH	SF181614B

Technical Specifications

Overall Dimensions (L x W x H)	in	43.7" x 21.3" x 36.4"
Number, Size of Tires	Front, Rear	2, 8" (2, 8")
Suggested User Weight ¹	lb	Max. 250
Weight of Heaviest Part ²	lb	33
Weight ³	lb	110
Each Battery Module Weight	lb	10
Battery Size	Type x PCS	12V14AH x 2
Charger	Amp	Off board (2A)
Output Power of Motor	HP	0.45
Maximum Speed	MPH	4
Maximum Climbing Angle/ Safe Climbing Angle ⁴	%(°)	10(6)
Range (Per Charge) ⁵	Mile	8.8
Turning Radius	in	46"
Ground Clearance	in	1.5"
Handle Bar	Type	T-Bar
Model #		7A
Standard Seat	W x D x BH	SF181614B

Note: ¹ On level Road. ² Without basket and seating. ³ With basket, seat & batteries. ⁴ Base on user weight of 200LB. ⁵ Based on user weight of 165LB. ⁶ After the batteries and mechanical moving parts fully break in.



Sunrunner 3

The Sunrunner 3 is a mid-size three-wheel scooter with a 300 lb. weight capacity. The unit comes standard with a deluxe captain's seat that can accommodate most users. It has an adjustable tiller for added user comfort and a front head light for added convenience for evening operation. It is available optionally with a full lighting package adding turn indicators and hazard lights. The unit has a 3in. ground clearance for operation on multiple terrains.

- + 300lb Weight Capacity
- + Slide & Swivel Seating
- + Adjustable Tiller
- + Optional Rear Light Package
- + Adjustable Seat Height



● Burgundy ● Blue

Model Number: 888B-3
Medicare Billing Code: K0806

HCPCS codes provided should not be considered as legal advice, and does not guarantee reimbursement.



Sunrunner 4

The Sunrunner 4 is a mid-size four-wheel scooter with a 300 lb. weight capacity. The unit comes standard with a deluxe captain's seat that can accommodate most users. It has an adjustable tiller for added user comfort and a front head light for added convenience for evening operation. It is available optionally with a full lighting package adding turn indicators and hazard lights. The unit has a 3in. ground clearance for operation on multiple terrains.

- + 300lb Weight Capacity
- + Slide & Swivel Seating
- + Adjustable Tiller
- + Optional Rear Light Package
- + Adjustable Seat Height



● Burgundy ● Blue

Model Number: 888B-4
Medicare Billing Code: K0899

HCPCS codes provided should not be considered as legal advice, and does not guarantee reimbursement.

Technical Specifications

Overall Dimensions (L x W x H)	in	44" x 23" x 42"
Number, Size of Tires	Front, Rear	1, 9" (2, 10")
Suggested User Weight ¹	lb	Max. 300
Weight of Heaviest Part ²	lb	-
Weight ³	lb	165
Each Battery Module Weight	lb	24
Battery Size	Type x PCS	12V35AH x 2
Charger	Amp	Off board (5A)
Output Power of Motor	HP	1
Maximum Speed	MPH	5
Maximum Climbing Angle/ Safe Climbing Angle ⁴	%(°)	14(8)
Range (Per Charge) ⁵	Mile	25
Turning Radius	in	37"
Ground Clearance	in	3"
Handle Bar	Type	T-Bar
Model #		888B-3
Standard Seat	W x D x BH	SV191821B

Technical Specifications

Overall Dimensions (L x W x H)	in	48" x 23" x 43"
Number, Size of Tires	Front, Rear	2, 10" (2, 10")
Suggested User Weight ¹	lb	Max. 300
Weight of Heaviest Part ²	lb	-
Weight ³	lb	180
Each Battery Module Weight	lb	24
Battery Size	Type x PCS	12V35AH x 2
Charger	Amp	Off board (5A)
Output Power of Motor	HP	1
Maximum Speed	MPH	5
Maximum Climbing Angle/ Safe Climbing Angle ⁴	%(°)	14(8)
Range (Per Charge) ⁵	Mile	25
Turning Radius	in	50"
Ground Clearance	in	3"
Handle Bar	Type	T-Bar
Model #		888B-4
Standard Seat	W x D x BH	SV191821B

Note: ¹ On level Road. ² Without basket and seating. ³ With basket, seat & batteries. ⁴ Base on user weight of 200LB. ⁵ Based on user weight of 165LB. ⁶ After the batteries and mechanical moving parts fully break in.



Sprinter XL4

The Sprinter XL4 is a heavy duty four-wheel scooter with a 350 lb. weight capacity. The unit comes standard with a deluxe captain's seat that can accommodate most users. It has an adjustable tiller for added user comfort and a front head light for added convenience for evening operation. It is available optionally with a full lighting package adding turn indicators and hazard lights. The 12in. solid wheels give the unit a 3in. ground clearance for operation on multiple terrains.

- + 350lb Weight Capacity
- + Slide & Swivel Seating
- + Adjustable Tiller
- + Optional Rear Light Package
- + Adjustable Seat Height



● Burgundy

Model Number: 889B-4
Medicare Billing Code: K0899

HCPCS codes provided should not be considered as legal advice, and does not guarantee reimbursement.

Technical Specifications

Overall Dimensions (L x W x H)	in	52" x 23" x 30-42"
Number, Size of Tires	Front, Rear	2, 12" (2, 12")
Suggested User Weight ¹	lb	Max. 350
Weight of Heaviest Part ²	lb	108
Weight ³	lb	181
Each Battery Module Weight	lb	32
Battery Size	Type x PCS	12V33-40AH x 2
Charger	Amp	Off board (5A)
Output Power of Motor	HP	1
Maximum Speed	MPH	5
Maximum Climbing Angle/ Safe Climbing Angle ⁴	%(°)	14(8)
Range (Per Charge) ⁵	Mile	25
Turning Radius	in	54"
Ground Clearance	in	3"
Handle Bar	Type	T-Bar
Model #		889B-4
Standard Seat	W x D x BH	SV191821B

Note: ¹ On level Road. ² Without basket and seating. ³ With basket, seat & batteries. ⁴ Base on user weight of 200LB. ⁵ Based on user weight of 165LB. ⁶ After the batteries and mechanical moving parts fully break in.



Enduro XL3⁺

The Enduro XL3⁺ is a heavy duty three-wheel scooter with a 500 lb. weight capacity. It has a delta style tiller which houses the controls for the lighting system, turn signals, wig wag and horn. It has a dual purpose wig wag which allows the user to push forward or pull back on the drive lever depending on dexterity for ease of driving. The unit comes standard with a 20in. seat that can accommodate most riders, and has swivel and recline features for added comfort and ease of use.

- + 500lb Weight Capacity
- + Slide & Swivel Seating
- + Full Suspension System
- + Full Lighting Package
- + Adjustable Tiller



● Burgundy

Model Number: 778XLSBN
Medicare Billing Code: K0899

HCPCS codes provided should not be considered as legal advice, and does not guarantee reimbursement.

Technical Specifications

Overall Dimensions (L x W x H)	in	53" x 28" x 52"
Number, Size of Tires	Front, Rear	1, 12" (2, 12")
Suggested User Weight ¹	lb	Max. 500
Weight of Heaviest Part ²	lb	104
Weight ³	lb	222
Each Battery Module Weight	lb	50
Battery Size	Type x PCS	12V50-75AH x 2
Charger	Amp	Off board (8A)
Output Power of Motor	HP	2
Maximum Speed	MPH	5
Maximum Climbing Angle/ Safe Climbing Angle ⁴	%(°)	14(8)
Range (Per Charge) ⁵	Mile	25-39
Turning Radius	in	49"
Ground Clearance	in	4"
Handle Bar	Type	Delta
Model #		778XLSBN
Standard Seat	W x D x BH	SV201820B



Enduro XL4⁺

The Enduro XL4⁺ is a heavy duty four-wheel scooter with a 500 lb. weight capacity. It has a delta style tiller which houses the controls for the lighting system, turn signals, wig wag and horn. It has a dual purpose wig wag which allows the user to push forward or pull back on the drive lever depending on dexterity for ease of driving. The unit comes standard with a 20in. seat that can accommodate most riders, and has swivel and recline features for added comfort and ease of use. It also features a full suspension to allow for smoother ride operation on most terrains.

- + 500lb Weight Capacity
- + Slide & Swivel Seating
- + Full Suspension System
- + Full Lighting Package
- + Adjustable Tiller

● Burgundy

Model Number: 889XLSBN
Medicare Billing Code: K0899

HCPCS codes provided should not be considered as legal advice, and does not guarantee reimbursement.

Technical Specifications

Overall Dimensions (L x W x H)	in	56" x 28" x 51"
Number, Size of Tires	Front, Rear	2, 12" (2, 12")
Suggested User Weight ¹	lb	Max. 500
Weight of Heaviest Part ²	lb	115
Weight ³	lb	240
Each Battery Module Weight	lb	50
Battery Size	Type x PCS	12V50-75AH x 2
Charger	Amp	Off board (8A)
Output Power of Motor	HP	2
Maximum Speed	MPH	5
Maximum Climbing Angle/ Safe Climbing Angle ⁴	%(°)	14(8)
Range (Per Charge) ⁵	Mile	25-39
Turning Radius	in	72"
Ground Clearance	in	4"
Handle Bar	Type	Delta
Model #		889XLSBN
Standard Seat	W x D x BH	SV201820B

Note: ¹ On level Road. ² Without basket and seating. ³ With basket, seat & batteries. ⁴ Base on user weight of 200LB. ⁵ Based on user weight of 165LB. ⁶ After the batteries and mechanical moving parts fully break in.

Accessories



USB Charger

Keep your devices charged with our USB charging adapter. Safely charge a phone or tablet while on the go from your scooter's charging port.

P/N: 206226-00100



Mesh Oxygen Tank Holder

A flexible way to store your oxygen tank while on the go. The mesh bag is light weight for easy installation on your product.

P/N: OXY-03



Elevating Leg Rests

An adjustable option to make your power Chair experience more comfortable. Our elevating leg rest's adjust to many angles and can swing away for easy entrance or exit.

P/N: ELR-888WA



Safety Flag

Be seen while using your product. The safety flag is 4.5 feet tall to make sure you can be seen while traveling around town.

P/N: SAF-9002



Metal Oxygen Tank Holder

A flexible way to store your oxygen tank while on the go. The mesh bag is light weight for easy installation on your product.

P/N:0210 (10inch length)
P/N:0216 (16inch length)



Swing Away Arm Rest

Our swing away joystick mount allows for extra space when moving your power Chair close to a table or desk.

P/N: 1343-00709-01-A3

Scooter End-to-End

All specifications are subject to change without notice. Shoprider Mobility Products, Inc. Reserves the right to any changes to the unit including A.) Recommended weight of user B.) Recommended ground surface and terrain use C.) Battery capacity and conditions D.) Type of charger E.) Ambient Temperature recommendations F.) Recommended method of driving G.) If the battery and mechanical moving parts are fully broken in.



Description	Units	Echo Folding	Echo 3	Hero	Scootie	Dasher 3	Escape	Sunrunner 3	Sunrunner 4	Sprinter XL4	Enduro XL3 ⁺	Enduro XL4 ⁺
		Folding	One Piece Frame	One Piece Frame	Length Adjustable	Take Apart	Take Apart	Mid-Size	Mid-Size	Heavy Duty	Heavy Duty	Heavy Duty
Overall Dimensions (L x W x H)	in	39.5" x 17.5" x 35"	37" x 21" x 36"	39" x 22" x 35"	37-40" x 17" x 34"	41.5" x 21.5" x 35"	43.7" x 21.3" x 36.4"	44" x 23" x 42"	48" x 23" x 43"	52" x 23" x 30-42"	53" x 28" x 52"	56" x 28" x 51"
Number, Size of Tires	Front, Rear	2, 6" (2, 8")	1, 8" (2, 8")	1, 8" (2, 8")	2, 7" (2, 8")	1, 8" (2, 8")	2, 8" (2, 8")	1, 9" (2, 10")	2, 10" (2, 10")	2, 12" (2, 12")	1, 12" (2, 12")	2, 12" (2, 12")
Suggested User Weight ¹	lb	Max. 200	Max. 250	Max. 250	Max. 250	Max. 250	Max. 250	Max. 300	Max. 300	Max. 350	Max. 500	Max. 500
Weight of Heaviest Part ²	lb	59	35	40	25	25	33	N/A	N/A	108	104	115
Weight ³	lb	59	70	73	82	87	110	165	180	181	222	240
Each Battery Module Weight	lb	7	15	15	24	22	10	24	24	32	50	50
Battery Module Capacity	Volt_Amps x Pcs	12V10Ah x 2	12V10Ah x 2	12V10Ah x 2	12V14Ah x 2	12V14Ah x 2	12V14Ah x 2	12V35Ah x 2	12V35Ah x 2	12V33-40Ah x 2	12V50-75Ah x 2	12V50-75Ah x 2
Charger		Off Board (2A)	Off Board (2A)	Off Board (2A)	Off Board (2A)	Off Board (2A)	Off Board (2A)	Off Board (5A)	Off Board (5A)	Off Board (5A)	Off Board (8A)	Off Board (8A)
Output Power of Motor	HP	.45	.45	.45	.75	.5	.45	1	1	1	2	2
Maximum Speed ²	MPH	3.75	3.8	3.8	4	3.8	4	5	5	5	5	5
Maximum Climbing Angle/ Safe Climbing Angle	%/(°)	10(6) ⁵	10(6) ⁴	10(6) ⁴	10(6) ⁴	10(6) ⁴	10(6) ⁴	14(8) ⁴	14(8) ⁴	14(8) ⁴	14(8) ⁴	14(8) ⁴
Range Per Charge ⁶	Mile	8.7 - 10	8	8	10	7.5	8.8	25	25	25	25-39	25-39
Turning Radius	in	36"	32"	33"	32"	33"	46"	37"	50"	54"	54"	72"
Ground Clearance	in	1.5"	1.5"	1.5"	1.5"	1.5"	1.5"	3"	3"	3"	4"	4"
Handle Bar	Type	T-Bar	T-Bar	T-Bar	T-Bar	T-Bar	T-Bar	T-Bar	T-Bar	T-Bar	Delta	Delta
Model Number	P/N	FS777	SL73	SL73N	TE-787NA	GK83	7A	888B-3	888B-4	889B-4	778XLSBN	889XLSBN

Note:
¹ On level Road. ² Without basket and seating. ³ With basket, seat & batteries. ⁴ Base on user weight of 200LB. ⁵ Based on user weight of 165LB. ⁶ After the batteries and mechanical moving parts fully break in.


SHOPRIDER
A TRUSTED NAME IN MOBILITY


21184 S. Figuerod St
Carson, CA. 90745
800.743.0772



 [shopridermobility](#)

 [shopridermobility](#)

 [shopridermobile](#)

 [shoprider.com](#)