



For easy ordering and information, please contact Customer Care Toll Free at 1(800) 743-0772 Fax: 1(877) 797-7081

6Runner 10

Name: _____
 Phone Number: _____
 Acct. Number: _____
 PO Number: _____
 Address: _____
 City: _____ State: _____ Zip Code: _____

Ship to: _____
 Phone Number: _____
 Address: _____
 City: _____ State: _____ Zip Code: _____

Base Scooter	Description	MSRP
<input type="checkbox"/> 888WNLM-RED	6Runner 10 Red	K0823 \$5,595.00
	6Runner 10 w/ Pan Seat Upgrade	K0822
BATTERY		
BAT-U1	Battery U1	x2 Incl.
SEATS (W"xD"xH")		
SV191821B	19X18 Captain Seat BLK w/o Recline	Incl.
ARM REST ASSY		
104300-00106-A3	(L) Flipback Height Adj. Armrest	Incl.
104400-00106-A3	(R) Flipback Height Adj. Armrest	Incl.
DRIVE TIRES		
?	Black 10" Flat Free Drive Tires	x2 Incl.
CASTORS		
1351-00800-01	7" Front and Rear Castor Assy	x4 Incl.
?	7" Black castor tire	Incl.
FOOT REST		
2329-00100-11	Flip Up Foot Plate	Incl.
CONTROLLER		
8301-00200-06	Shark Joystick Controller	Incl.
1361-00204-04	Shark Main Control Box	Incl.
CHARGER		
109201-88124	5 Amp Charger Off-Board	Incl.
ACCESSORIES		
300909-03	Safety Belt	Incl.
Total		_____

All specifications are subject to change without notice. Shoprider Mobility Products, Inc. reserves the rights to any changes to the unit.
 a) Weight of user b) Ground surface and terrain c) Battery capacity and conditions d) Type of charger e) Ambient temperature f) The method of driving g) If the battery and mechanical moving parts are fully broken in.

SPECIFICATIONS

OVERALL DIMENSIONS

Base Length 42"
 Base Width 23"
 Height w/o headrest..... 46"

Ground to Seat Top

Ground Clearance 2"

WEIGHTS

User Weight Capacity 300 lbs
 Total Unit Weight 191 lbs

Battery Weight Module 23 lbs

Tires

Number, Size of Front Tire 2, 7"
 Number, Size of Mid Wheel 2, 10"
 Number, Size of Rear Tire 2, 7"

Batteries

BAT - U1 Module Capacity 12v32Ahx2

Charger

Charger-On Board 5 amp

Motor Power Output

Horsepower60 hp x2

Driving Specs

Turning Radius 15.75"
 Range 10 miles
 Maximum Speed 5 mph
 Maximum Climbing Angle 10% / 6°



Replacement Parts & Accessories

Part Number	Description		MSRP
BATTERY			
BAT-U1	Battery U1	x2	\$140.00ea
SEATS (W"xD"xH")			
<input type="checkbox"/> SV191821B	19X18 Captain Seat BLK w/o Recline		\$280.00
<input type="checkbox"/> SV201820RB	20x18x20 Recline Blk Pillow Top		\$370.00
<input type="checkbox"/> SV2121211/2RG	21x21x211/2 Recline Capt Seat-Grey Vinyl		\$400.00
<input type="checkbox"/> SV2121211/2RB	21x21x211/2 Recline Capt Seat-Blk Vinyl		\$400.00
<input type="checkbox"/> SV191820RG	19x18 Recline High-Back Grey Vinyl		\$400.00
<input type="checkbox"/> SP181821RG	18x18 Pan seat w/Capt. Back Recline		\$400.00
ARM REST ASSY			
<input type="checkbox"/> 104300-00106-A3	(L) Flipback Height Adj. Armrest	x2	\$140.00
<input type="checkbox"/> 104400-00106-A3	(R) Flipback Height Adj. Armrest	x2	\$140.00
<input type="checkbox"/> 104310-00102	Arm pad w/ height-adjustment tube	x2	\$50.00ea
DRIVE TIRES			
<input type="checkbox"/> TSR26085 Rear	Grey 10" Flat Free Drive Tires	x4	\$50.00ea
<input type="checkbox"/> No Part # Yet	Blk 10" Flat Free Drive Tires	x4	\$50.00ea
CASTORS			
<input type="checkbox"/> 1351-00800-01	7" Front & Rear Castor Assy		\$140.00ea
<input type="checkbox"/> 105101-69000	7" castor tire only		\$60.00ea
<input type="checkbox"/> No Part # Yet	7" castor tire only Blk		\$60.00ea
FOOT REST			
<input type="checkbox"/> 2329-00100-11	Flip Up Foot Plate		\$150.00
<input type="checkbox"/> ELRJET-R(K)	Elevating Leg Rest Kit-Red		\$232.00
<input type="checkbox"/> ELRJET-B(K)	Elevating Leg Rest Kit-Blk		\$232.00
<input type="checkbox"/> 202934-00104	(L) Leg Rest Hanger-ELR/SAFR		\$32.00
<input type="checkbox"/> 202934-00105	(R) Leg Rest Hanger-ELR/SAFR		\$32.00
<input type="checkbox"/> 800907-007-00	(L)Swing Away Foot Rest		\$80.00
<input type="checkbox"/> 800907-007-01	(R)Swing Away Foot Rest		\$80.00
CONTROLLER			
<input type="checkbox"/> 8301-00200-06	Shark Joystick Controller		\$360.00
<input type="checkbox"/> 1361-00204-04	Shark Main Control Box		\$776.00
CHARGER			
<input type="checkbox"/> 109201-88124	5 Amp Charger Off-Board		\$230.00
ACCESSORIES			
<input type="checkbox"/> O210(K)	Oxygen Tank Holder (Metal)		\$150.00
<input type="checkbox"/> WH-0185(K)	Walker Holder		\$150.00
<input type="checkbox"/> OXY-03(K)	Mesh Oxygen Tank Holder		\$60.00
<input type="checkbox"/> RB-0125(K)	Rear Basket		\$170.00
<input type="checkbox"/> 300909-03	Safety Belt		\$34.00

Total _____