




EZ BATHING ACCESSORIES

Bathe Easy!

With sink, bath, and shower accessories from EZ-ACCESS®. Our complete line of quality bathing accessories can help make personal hygiene more convenient, safe, and enjoyable in the comfort of your own home.



Enjoy bathing the way it's meant to be – soaking wet! The **EZ-BATHE® Body Washing Basin** allows for a refreshing in-bed bath. The tub is unrolled and inflated around the patient, allowing for a convenient and safe bathing experience for both the patient and caregiver.

- Fits standard and hospital beds
- Accommodates individuals up to 6' 2"
- Features a double-tube design to help prevent splashes and spills
- Inside dimensions of inflated tub are 71" long, 22" wide, and 13½" deep
- Made of heavy-duty vinyl
- Latex free
- One-year limited warranty

Includes:

- Wet & dry vacuum to inflate/deflate the tub
- 25' hose with hand-held shower head attaches to sink for water supply
- 25' drain hose with on/off valve



Our exceptionally easy-to-use **EZ-SHAMPOO® Hair Washing Tray** is contoured to fit the neck comfortably and rest on the client's shoulders for convenient shampooing in a seated position.

- Compact and lightweight design
- Reduces awkward lifting and bending for caregiver
- Raised edge prevents water from spilling over
- 35" stabilizing strap included
- Designed for use with chair or wheelchair
- Made of durable plastic and designed for years of trouble-free use
- One-year limited warranty

The **EZ-SHAMPOO® Hair Washing Basin** allows for a refreshing, in-bed shampoo. Comfortable and convenient, the inflatable basin provides secure, form fitting support that cushions the head and shoulders and includes an adjustable, built-in headrest. A caregiver can quickly shampoo and rinse a patient's hair without awkward lifting and bending. An attached one-way drain hose easily and conveniently removes suds and water.

- Inflatable and lightweight for portability and convenient storage
- Features a double-tube design to help prevent splashes and spills
- Flexible drain hose with custom stopper
- 24" wide x 20" long x 8" deep, inflated size
- 6½" neck opening (one size fits most)
- Made of heavy-duty vinyl
- Latex free
- One-year limited warranty



Designed to be used where a water source is unavailable, the **EZ-SHOWER** is an overhead bedside shower and is the ideal accessory to the **EZ-SHAMPOO® Hair Washing Basin**.

- Holds up to 2½ gallons of warm tap water – ample enough for a complete shampoo and rinse
- Flexible hose directs a gentle stream of water exactly where it's needed
- Compact and easy to store when not in use
- Includes a flexible 2' hose with an on/off switch to control water flow
- Built-in handle to carry to the bedside
- Hangs on a bedpost or an I.V. pole
- Made of heavy-duty vinyl
- Latex free
- One-year limited warranty

customerservice@ezaccess.com | 1.800.451.1903 | www.ezaccess.com