



# Opus Rollator

**300 lbs**  
WEIGHT CAPACITY



ITEM #	DESCRIPTION
960GG	Euro Style 4 Wheel Carbon Fiber Rollator (Graphite Grey)

## FEATURES

- 1 Lightweight durable carbon fiber frame for easy transport.
- 2 Backrest is height and depth adjustable for optimal comfort.
- 3 Cross-brace frame for added strength and safety.
- 4 Quick and easy to fold up.
- 5 Large removable tote bag for convenient storage included.
- 6 Large 10" Chrome Front Wheels for easy maneuvering for indoor/outdoor use.
- 7 Clamp lock holds rollator closed when it's folded.

## SPECIFICATIONS

Weight	14.5 lbs
Handle Height	32.5" - 36.25"
Height Adjustment	Handles Only
Seat to Floor Height	21"
Adjustable Backrest	4" Height Adjustment and 1.5" Depth Adjustment
Seat Dimensions	18" x 8.75"
Unit Dimensions	24.5" W x 29" L
Frame	Carbon Fiber
Weight Capacity	300 lbs

## WARRANTY

Lifestyle Mobility Aids warrants this product to be free of defects in materials and workmanship for the lifetime of the original consumer. We offer this **Limited Lifetime Warranty** to affirm our confidence that the materials and workmanship of our products will result in years of dependable service.

This warranty is conditioned upon the proper use of this product in accordance with the instructions and manuals, including proper use and maintenance. This warranty does not cover non-durable components, such as rubber accessories, which are subject to normal wear and will require periodic replacement. In the event warranty service is required, we will determine whether the product is to be repaired or replaced.

For questions or to obtain warranty service, please contact an authorized Lifestyle Mobility Aids dealer.