



Adjustable hold-down arm safely secures wheelchair for transport.

Simply Ideal for Wheelchair Transport

Lightweight, powered lift transports a manual folding wheelchair easily. Features a 100 lbs. capacity with simple, one-switch operation. Adjustable hold-down arm safely secures the chair. Lift folds when not being used. Swing-Away option available.

Features & Benefits

- One-switch operation
- Adjustable hold-down arm
- Power lift, manual fold
- No chair modifications required
- Folds when not in use
- Ball mount standard; upgrades available
- Manual crank backup and license plate mount included
- Control box with license plate mount includes easy switch, key lock, and courtesy light
- Swing-Away option available

Specifications	AL003 Tilt-N-Tote
Accommodates	Manual folding wheelchairs
Power	Power lifting
Lifting Capacity	100 lbs. (45 kg)
Installed Weight	49 lbs. (22 kg)
Securement	Adjustable manual hold-down arm
Hitch	Class I or greater
Standard Hitch Height	5" – 26"
Included Accessories	Ball mount
Available Options	Battery pack, manual wheelchair covers
Warranty	3-year transferable



Optional AL105 Swing-Away allows easy access to vehicle's cargo area.

Optional AL205 Battery Pack avoids hard-wiring into vehicle or serves as emergency backup.

