

AL500 / AL500HD Universal Power Chair Lifts



Fits Virtually All Power Chairs

Harmar's AL500 and AL500 Heavy Duty accommodates virtually any power chair. Available in traditional 350 lbs. or heavy-duty 400 lbs. models.

Features & Benefits

- Easy, one-switch operation, license plate light, key switch and courtesy light
- (4) self-tensioning Q'Straint retractors; no chair modifications required
- Manual crank backup and license plate mount included
- Swing-Away option available



Optional AL105 Swing-Away Shown

Swing-Away Options

This easy to operate option allows an unoccupied outside lift to be easily pivoted away for clear access to the trunk, lift gate or tailgate.

- **AL105 Swing Away**- 400 lbs capacity, 38 lbs install weight
- **AL105L Light Weight Swing-Away**-350 lbs capacity, 25 lbs install weight. Lighterweight offers greater handling/ steering and fuel consumption.

| Specifications | AL500 | AL500 HD Heavy Duty |
|------------------|---|---|
| Accommodates | Power wheelchairs | Power wheelchairs |
| Platform | 27.5" x 39" aluminum deck | 27.5" x 39" aluminum deck |
| Lifting Capacity | 350 lbs. (159 kg) | 400 lbs (181 kg) |
| Installed Weight | 75 lbs. (34 kg) | 77 lbs. (35 kg) |
| Securement | (4) Self-tensioning Q'Straint retractors | (4) Self-tensioning Q'Straint retractors |
| Hitch | AL123 class III hitch adapter and AL107 hitch post; class II optional | AL123HD class III hitch adapter and AL107 hitch post included |
| Options | Battery pack, Swing-Away, Power wheelchair covers | Battery pack, Swing-Away, Power wheelchair covers |
| Warranty | 3-year transferable | 3-year transferable |



800-833-0478 | www.harmar.com

