

AL600 Hybrid Platform Lift



NEW AND IMPROVED DESIGN!
Wires and chain now covered

Drive on, drive off platform

Easy to operate remote control

Optional wireless remote

Lift connects to existing third row seating!


MADE IN USA

FEATURES

- Carries virtually all standard power chairs & scooters
- Weatherproof, remote hand control (wireless option available)
- No drilling required in most vehicles
- No scooter or power chair modifications
- Most compact unit on the market!
- Allows continued use of second row seating
- Manual crank backup
- Compatible with most minivans and SUVs

SPECIFICATIONS

Lifting capacity:	350 lbs
Platform size:	28.5" x 38"-50"
Headroom required:	32"
Depth required:	43" maximum
Platform height:	5.5" maximum (dead space)
Vertical reach:	33.5"
Warranty:	3-year transferable

DESCRIPTION

Harmar's AL600 combines the simplicity of a drive-on, drive-off platform with all the benefits of carrying your mobility device inside a vehicle.

With only 43" of depth and 32" of headroom required, this is **the most compact hybrid lift on the market today**. Though compact, it still manages to offer the largest travel distance in its class, making this lift compatible with most enclosed vehicles, including SUVs! Installation utilizes existing seat mounts, requiring no drilling or further vehicle modifications. To top it all off, the incredibly compact design allows the second row seating to remain in place!