

DIMENSIONS		FDA, TÜV, CE Approved
Trolley mode:	H 28"	71 cm
	L 15.3"	38.9 cm
	W 16.5"	41.9 cm
Drive mode:	H 35"	89 cm
	L 47"	120 cm
	W 22"	56 cm
WEIGHT		
Front half:	26.4 lbs	12 kg
Rear half:	37.5 lbs	17 kg
Assembled:	64 lbs	29 kg (w/o battery)
Battery:	4.1 lbs	1.85 kg
WEIGHT CAPACITY		
Maximum user weight:	264 lbs	120 kg
BATTERY		
48 Volts Lithium-ion, more than 20 km / 12 miles of continuous driving, only 4 hours on average for full charge, charge on-board and off-board		
SPEED		
Most countries, 4 mph/ 6.4 km per hour*		
* According to the local regulation		



Separates into
2 parts in seconds



Easy to wheel
as a trolley



Compact to
full size in 10 sec

ATTO

The best folding scooter in the world