

Stealth Products Cushion Specifications

| Cushion Name | Weight Capacity | Weight of Cushion | Cushion Height | Width Range | Depth Range | HCPCS Code | MSRP |
|-------------------|-----------------|-------------------|----------------|-------------|-------------|----------------|--|
| Simplicity G | 350 lbs | 2.2 lbs | 2" | 16" - 24" | 16" - 20" | E2601 E2602 | \$110.00 \$165.00 |
| Solution SPP | 350 lbs | 5.4 lbs | 3" | 10" - 24" | 10" - 24" | E2607 E2608 | \$375.00 \$430.00 |
| Spectrum Foam SPP | 650 lbs | 3 lbs | 3" | 16" - 32" | 16" - 28" | E2607 E2608 | \$460.00 \$530.00 \$590.00 >24"w |
| Spectrum Gel SPP | 350 lbs | 3.5 lbs | 3" | 14" - 24" | 16" - 24" | E2607 E2608 | \$460.00 \$530.00 |
| Tru Comfort 2 SPP | 450 lbs | 4 lbs | 4.5" | 16" - 22" | 14" - 24" | E2607 E2608 | \$400.00 \$500.00 |

Disclaimer: HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage, and bundling guidelines may apply. Stealth Products, LLC's review of Medicare guidance suggests that accessories acquired at zero-dollar net cost may not be eligible for separate Medicare reimbursement despite being separately coded and eligible to be unbundled at the time of initial purchase. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Stealth Products, LLC recommends that providers seek competent legal advice to assure a compliant Medicare billing protocol.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.



Solutions For Life

stealthproducts.com



Solutions For Life



POWERED BY coolcore®

stealthproducts.com

Stealth Products Cushions Color Guide

To make your cushion selection process a little easier, Stealth Products Cushions will all be marked with a color coded pvc patch. The color of the patch will identify the function and HCPC code of the cushion.



E2601/E2602



E2603/E2604



E2605/E2606



E2607/E2608



E2622/E2623



E2624/E2625



E2609

SIMPLICITY G



E2601/E2602

Features:

- Contoured, high-density molded foam base
- Coccyx and sacral relief
- Moisture resistant and easy to clean Stealth-Tek incontinence cover.
- Sizes available - 16" - 24"W x 16" - 20"D
- Mesh cover optional

SPECTRUM FOAM SPP



E2607/E2608

Features:

- Lateral and pre-ischial shelf contoured, high-density molded foam base
- Multi-layer, viscoelastic foam insert under at risk bony areas
- CoolCore® technology cover standard.
- Urethane liner for incontinence, 4-Way Stretch Poly liner optional
- Sizes available - 16" - 32"W x 16" - 28"D
- Custom sizes available
- Optional Stealth-Tek, Mesh

SOLUTION SPP



E2607/E2608

Features:

- Lateral and pre-ischial shelf contoured, high-density molded foam base
- Twin-cell, aqueous gel insert under at risk bony areas
- CoolCore® technology cover standard.
- Urethane liner for incontinence
- Sizes available - 10" - 24"W x 10" - 24"D
- Custom sizes available
- Optional Stealth-Tek or Mesh cover.

SPECTRUM GEL SPP



E2607/E2608

Features:

- Lateral and pre-ischial shelf contoured, high-density molded foam base
- Twin-cell, lightweight fluid insert under at risk bony areas
- CoolCore® technology cover standard.
- Urethane liner for incontinence 4-Way Stretch Poly liner optional
- Sizes available - 14" - 24"W x 16" - 24"D
- Custom sizes available
- Optional Stealth-Tek, Mesh

TRU-COMFORT 2 SPP

E2607/E2608



Convenient side pouch on both sides of the cushion with easy pull rubber loop zipper handles. Pouch is large enough to fit a larger format mobile phone.


Features:

- Deep contoured, high-density molded foam base
- CoolCore Technology cover standard.
- Easy pull rubber loop zipper handles
- Convenient side pouches on both sides of the cushion
- Sizes available - 16" - 22"W x 14" - 24"D
- Custom sizes available
- 4-Way Stretch Poly liner optional

POWERED BY coolcore®

What is CoolCore®

- Chemical Free material
- UPF/Flammability rated
- 4-Way Stretch nylon interlock with Elastane
- Patented combination of chemical-free blends that deliver 3 distinct functions: **Moisture Wicking** to **Moisture Transportation** to **Moisture Evaporation**.
- **CoolCore® Technology is exclusive to Stealth Products.**



Year Round Temperature Control

CoolCore fabrics wick moisture away from the body in every direction keeping you dry and comfortable no matter what season it is.

Breakthrough Technology

Research driven, patented fabrics from Coolcore are chemical-free and yield remarkable benefits. They cool better, last much longer and offer protection from the sun's harmful rays (UPF).