

## **EW-BUGEYE Specifications**

COLORS: RED or LIGHT BLUE

MAXIMUM SPEED: 15 MPH (varies based on riders weight and driving terrain)

TRAVEL RANGE: 40 MILES (varies based on riders weight and driving terrain)

Bluetooth: Yes, Advanced Bluetooth system with sound

MOTOR: 500 watt Trans-Axel Motor

WEIGHT CAPACITY: 350 LBS

BATTERIES: 48 Volt, 20 AH

REAR VIEW MIRRORS: Yes, Two Mirrors

BRAKES: Front and Rear Drum Brakes

BRAKE LEVERS: Front and Rear

STORAGE: Rear Box Storage & Under Seat Storage both secure lockable

VARIABLE TWIST THROTTLE: Yes

SEAT SIZE: 14 "Wide and 23" Deep

SEAT BASE HEIGHT: From Floor Height 30"

SMART CHARGER: Yes

BAG HOOK: Yes

PARKING BRAKE: Yes

Front Storage Area: Two Storage Bins

**REVERSE SWITCH: Yes** 

**ELECTRIC HORN: Yes** 

GUAGES: Yes, Speedometer with Battery Gauge

3 Speed Settings Low, Medium and High

CLIMBING: 12 degrees

LIGHTING PACKAGE: Turn Signals, Brake Lights, High/Low Beam Headlight, Rear Light

KEYLESS KEY FOB: Yes, Power on/off

ANTI-THEFT ALARM SYSTEM: Yes, Keyless Remote

BATTERY INDICATOR GAUGE: Yes

TIRE SIZE: Front and Rear Tire Size 15"/3"

REAR REFLECTORS: Yes

ALUMINUM RIMS: Yes

SCOOTER LENGTH: 70"

SCOOTER WIDTH: 28"

SCOOTER HEIGHT: 40"

DECK HEIGHT: 12"

GROUND CLEARANCE: 6.5"

WHEEL BASE: 52"

FRONT SUSPENSION FORK: Yes

REAR SUSPENSION: Two Heavy Duty Shocks

SCOOTER WEIGHT: 209 LBS

SCOOTER SHIP WEIGHT: 335 LBS

WARRANTY 3 YEAR LIMITED