

Compare Luggie Models

Features	<u>Classic</u>	<u>Standard</u>	<u>Elite</u>	<u>Super</u>
Weight Capacity	250 lbs	250 lbs	320 lbs	360 lbs
Range*	10 Miles	10 Miles	15 Miles	11 Miles
Max Speed	4 mph	4 mph	4 mph	4.5 mph
Turning Radius	35.5"	35.5"	35.5"	35.5"
Ground Clearance	2.5"	2.5"	2.5"	2.5"
Battery	12AH	8.5AH	10.5AH	10.5AH
Lithium Battery	No	Yes	Yes	Yes
Battery Weight	20 lbs	5 lbs	5 lbs	5 lbs
Charging Time	2 hours	90 mins	90 mins	90 mins
Battery Indicator	Meter	Meter	LED	LED
Seat Width	14.5"	14.5"	15.7"	17"
Armrests Standard	No	No	Yes	Yes
Seat Height Options (20" & 22")	Yes	Yes	Yes	Yes
Overall Length: (unfolded)	39"	39"	39"	39"
Overall Length (folded)	26"	26"	28"	28"
Overall Width	18"	18"	18"	20.7"
Speed Control Dial	Yes	Yes	Yes	Yes
Extendable Handlebars	No	No	Yes	Yes
Adjustable Steering Heights	Yes	Yes	Yes	Yes
Auto Retractable Bumper	No	No	Yes	Yes
Total Weight (without battery)	51.7 lbs	51.7 lbs	51.7 lbs	67 lbs
Aluminum Alloy Frame	Yes	Yes	Yes	Yes
Rear Wheel Drive	Yes	Yes	Yes	Yes
24 Volt DC Magnetic Motor	Yes	Yes	Yes	Yes
Intelligent Braking System	Yes	Yes	Yes	Yes
Airline Safe	Yes	Yes	Yes	Yes