




LC-105 shown in  
Lexis Urethane® New Chestnut

Protect your investment with the Pride® Power Lift Recliner Care Kit



# Essential COLLECTION

The Essential Collection from Pride® Power Lift Recliners provides exceptional value and outstanding comfort. Several fabric options will match any décor while delivering rest and worry-free operation.

AMERICA'S #1  
POWER LIFT RECLINER

POWER  
RECLINERS

US: 1-800-800-8586 • CAN: 1-888-570-1113  
pridemobility.com |    

**PRIDE**  
MOBILITY  
LIVE YOUR BEST™

# essential COLLECTION

AMERICA'S #1  
POWER LIFT RECLINER

POWER  
LIFT  
RECLINERS



LC-250 shown  
in Cloud 9 Stone



	Medium	Medium
SPECIFICATIONS	LC-250	LC-105
	3-Position Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 5'2" - 5'10"
Weight Capacity	375 lbs.	325 lbs.
Seat-to-Floor	20"	18"
Seat Depth	20"	19"
Seat Width	22"	20"
Top of Back to Seat	26"	28"
Back Style	Blown Fiber Split Back	Sewn Pillow Back
Fully Padded Chaise	Yes	Yes
Heat & Massage	No	No
Seating Type	Total Comfort	Sinuous Spring
Distance From Wall	20"	17"
Overall Chair Width	34"	33"
Chair Weight	118 lbs.	96 lbs.
Warranty	2-Year Limited	1-Year Limited

## L-250 FABRIC SELECTION



### CLOUD 9



PACIFIC



STONE



WALNUT



BLACK CHERRY

OPTIONAL

### LEXIS STA-KLEEN®



CHESTNUT

## LC-105 FABRIC SELECTION



### MICRO-SUEDE



COCOA



SANDAL



SKY

OPTIONAL

### LEXIS URETHANE®



NEW CHESTNUT

## BENEFITS

- Easy-to-use hand control adjusts your chair to find your perfect position
- Large dual pockets offer ample storage space that is easily accessible
- Smooth and durable lift and scissor mechanism for long-lasting and quiet performance
- Removable back lets you adjust chair padding as needed
- Integrated battery backup ensures chair is operational when the power goes out<sup>1</sup>
- Energy-efficient power box equipped with self-testing electronics and 9' cord
- Lifetime warranty on steel lift frame
- Fire retardant (California specifications; bulletin 117-2013)

## OPTIONS\*

- Cleaning kit
- Non-skid leg levelers

\*Additional charges apply for all options.

<sup>1</sup>. Batteries not included.



Look for the Quick Ship Indicator for the shortest lead times available.



## 3-POSITION FULL RECLINE

The 3-position full recline is available on most Pride® Power Lift Recliner models and puts you in a position ideal for napping or relaxing.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

US: 1-800-800-8586 • CAN: 1-888-570-1113  
pridemobility.com |

**PRIDE**  
MOBILITY  
LIVE YOUR BEST™