



EW-36 Elite Spec Sheet

WEIGHT CAPACITY:

350 lbs.

POWER:

Electric

WATTS:

500 Watt

MOTOR TYPE:

Transaxle differential brushless motor

AMPS:

20 AH

VOLTS:

48 Volt

DIMENSIONS:

61" * 29" * 41"

BATTERIES:

48 volt lead acid maintenance free battery

TIRE SIZE: Front and Rear

15 x 3.5

CHARGER:

Smart charger included

SPEED:

Up to 18mph

Turning Radius:

50"

Ground clearance:

5"

DISTANCE:

Up to 45 miles per charge (distance varies by riders weight , terrain, road surface etc.)

HEADREST:

Yes adjustable headrest

THROTTLE TYPE:

Variable speed control- Twist Throttle

KEY START:

Yes

BRAKING SYSTEM:

Front and rear braking system

DRIVE SYSTEM:

Electromagnetic Brakes

HEAD LIGHT:

Yes

CARTON SIZE:

Shipped on a pallet

Weight Capacity:
350 lbs.

SHIPPING WEIGHT:
280 lbs

NET WEIGHT:
200 lbs

ELECTRIC LIGHT:
Yes

BATTERY INDICATOR:
Yes

STORAGE BASKET:
Yes

SHOCK ABSORBER:
Two rear shock absorber on the frame

ALARM:
Anti-theft alarm

WARRANTY:
Limited 3-Year Warranty