

Model EW-54 Spec Sheets

Scooter Weight Capacity 500Lbs.

Speed: Up to 15mph

Distance: 43 miles per charge (distance varies on terrain, riders weight and road surface etc.)

Front Seat: Type Adjustable chair with backrest

Front seat: Size 22 ¹/₂"x17 ¹/₂"x23" (LxWxH)

Front Seat: 27" From Ground

Rear Seat: Type No rear seat

Rear Seat: Size No rear seat

Steering Height 41" from Ground

Cargo Box Dimensions: 12"x20"x11" (LxWxH)

Power: Electric Watts: 700 Watt

Motor Type: Transaxle differential brushless Amps: 20 AH

Volts: 60 Volts

Dimensions: 72" x 27" x 65" (LxWxH)

Floor Clearance: 5"

Front Basket: No Front basket provided

Batteries: (5) Sealed Lead Acid Maintenance Free

Forward/Reverse switch: Yes

High/ Low Speed switch: (3) Speed switch

Front tire size: 13" Tire Type: 90/65-8

Rear tire size: 15" Tire Type: 120/50-9

Charger: 60-volt Smart Charger included

Throttle Type: Variable Speed Control – Twist throttle

Key Start: Yes

Braking System: Front and Rear drum brakes

Drive System: 700 Watt Brushless Motor

Headlight: 12V Bulb

Rear lights: 12V Bulb

Carton Size: Shipped Fully assembled on Pallet – 80"x36"x73" (LxWxH)

Scooter Weight: 333lbs. with batteries

Battery Indicator: Yes

Front: Front Shock Absorbers

Rear: Rear Shock Absorbers

Reflectors: Yes

Warranty: Limited – 3 Year Warranty