



EW-M31 COMPACT POWER CHAIR

Key Specifications:

Top Speed: 4 mph

Drive Range: 15.2 miles

Weight Capacity: 300 lbs.

Seat Type: 18.5" Swivel Captain Seat with Adjustable Headrest & Seat Height

Turning Radius: 31"

General Specifications:

Swing-Away joystick: Yes

Lap Belt: Yes

Variable Speed Dial: Yes

Optional Left or Right Hand Joystick: Yes

Brake: Electro-Mechanical

Self-Diagnostic Warning Lights: Yes

Drive Type: 2 Motor Rear-Wheel Drive

Grade Climbable: 8 degree

Footplate: Large flip-up Footplate

Adjustable Footplate: Yes, You can change the footplate vertically or horizontally

Ground Clearance: 2.4”
Adjustable seatback
Adjustable Armrest: Yes, width adjustable
Rear bag storage
Removable Seat: Yes
Seat Width: 18.5”
Seat Depth: 17”
Seat to Floor Height: 21”
Overall Width: 24”
Overall Length: 37.5”
Overall Height: 40”
Electric Horn: Yes
Battery Gauge: Yes
Free Wheel Levers: Yes
Joystick Options: Right or Left
Extendable Joystick Adjustment: Yes
Battery: (2) 12V 36AH
Battery Included: Yes
Battery Charger Type: Off Board 3A
Rear Wheels Size: 10”
Front Wheels Size: 7”
Solid Flat Free Wheels: Yes
Base Product Weight: 176.8 lbs.
Anti-Tip Wheels: Yes
Rear & Side Reflectors: Yes

