



The EW-M35 sets the standard for quality and portability. With its tubular design and reinforced frame structure, none other will compare. The sleek design, padded armrest provides the rider with both comfort and performance. The removable battery pack with wireless connection can be charged on or off the unit. Equipped with puncture proof tires, swivel seat and delta tiller are just a few features that makes EWheels Medical a world class leader in mobility.

Features include a maximum speed of 4 MPH, a 300 pound weight capacity, puncture proof tires, compact size for easy transportation and disassembles in seconds.

### EW-M35 Specifications

Capacity:  
285 lbs

Overall Length  
39.4"

Overall Width  
20"

Overall Height

35"

Weight of Heaviest Piece

33 lbs

Wheels:

Front 8"

Rear

8"

Weight w/ Batteries

88 lbs

Max Speed

4mph

Ground Clearance

1.2"

Grade Climbable

8 degree

Turning Radius

49"

Brake

Electromagnetic

Seat Type

Padded Foldable

Seat Width

17"

Motor Size

250w

Battery Size  
(2) 12V 12AH

Battery Weight  
20lbs

Travel Range  
9 miles

Electronics  
On/Off switch, Battery indicator