



EW-M41

Features include a maximum speed of 5mph, a weight capacity of 300 lbs., a range of 15 miles, a 22 AH battery pack for longer travel range, LED headlight, easy take apart mechanism, on board battery charger, and delta handlebar.

EW-M41 Specifications:

Overall Length
44"

Overall Width
23"

Overall Height
37"

Wheels:
Front 9"
Rear 9"

Weight w/ Batteries
121 lbs.

Max Speed
5mph

Weight Capacity
300 lbs.

Ground Clearance
2"
Grade Climbable
8 degree

Turning Radius
55"
Brake
Electromagnetic

Seat Type
Swivel Light Weight Foldable

Seat Width
18"

Motor Size
250 W

Battery Size
(2)12V 22AH

Battery Weight
34lbs.

Travel Range
15 Miles

Battery Charger
On Board

Electronics
On/Off Key Switch, Battery Indicator, Speed Control Knob