

EW-M47 HD Folding Power Wheelchair

Key Specs:

Seat Width: 21"

Top Speed: 5 mph

Battery Charge Distance: 12 miles

Weight Capacity: 350 lbs.

General Specifications:

Swing-Away joystick: Yes

Folded Size: 31"H x 32"D x 18" W

Lap Belt: Yes

Turning Radius: 33"

Optional Left or Right Hand Joystick: Yes

Brake: Electro-Mechanical

Frame Material: Aluminum

Top Speed: 5 mph

Adjustable Variable Speed: Yes

Self-Diagnostic Warning Lights: Yes

Heaviest Piece: 77 lbs.

Battery Charge Distance: 12 miles

Drive Type: 2 Motor Rear Wheel Drive

Grade Climbable: 8 degree
Foldable: Yes
Footplate: Large Footplate
Adjustable Footrest Height: Yes
Ground Clearance: 3.2"
Seat Width: 21"
Seat Depth: 15.6"
Seat to Floor Height: 21.75"
Overall Width: 27.4"
Overall Length: 41"
Overall Height: 40"
Electric Horn: Yes
Battery Gauge: Yes
Free Wheel Levers: Yes
Joystick Options: Right or Left
Extendable Joystick Adjustment: Yes
Frame Type: Aluminum
Battery: (2) 12V 36AH
Battery Included: Yes
Battery Charger Type: Off Board 5A
Rear Wheel Size: 13"
Front Wheel Size: 7.9"
Pneumatic Rear Wheel: Yes

Solid Front Wheel: Yes

Base Product Weight: 143 lbs.

Battery Pack Weight: 55 lbs.

Weight Without Batteries: 88 lbs.

Anti-Tip Wheels: Yes