## EW-M51 HEAVY DUTY POWER CHAIR

**General Specifications:** 

Swing-Away joystick: Yes Lap Belt: Yes Variable Speed Dial: Yes **Turning Radius: 31**" **Optional Left or Right Hand Joystick: Yes Brake: Electro-Mechanical** Top Speed: 5 mph Self-Diagnostic Warning Lights: Yes **Battery Charge Distance: 15 miles** Drive Type: 2 Motor Rear-Wheel Drive Grade Climbable: 8 degree Footplate: Large flip-up Footplate Adjustable Footplate: Yes, You can change the foot plate angle and height Ground Clearance: 3.4" Seat Type: 22" Swivel Captain Seat with Adjustable Headrest & Seat Height Adjustable Armrest: Yes, width adjustable **Removable Seat: Yes** 

**Removable Batteries: Yes** Seat Width: 22" Seat Depth: 21" Seat to Floor Height: 22" **Overall Width: 25**" **Overall Length: 44**" **Overall Height: 46**" **Electric Horn: Yes Battery Gauge: Yes** Free Wheel Levers: Yes Joystick Options: Right or Left Extendable Joystick Adjustment: Yes Battery: (2) 12V 50AH **Battery Included: Yes** Battery Charger Type: Off Board 5A Rear Wheels Size: 14" Front Wheels Size: 9" Solid Flat Free Wheels: Yes Base Product Weight: 224 lbs. Battery Pack Weight: 67 lbs. Weight Without Batteries: 157 lbs. Anti-Tip Wheels: Yes Rear & Side Reflectors: Yes

Weight Capacity: 400 lbs.