



3 wheel GB117Z 4 wheel GB147Z The Buzzaround XLSHD features compact portability with a heavy-duty weight capacity! The Buzzaround XLSHD features a robust 325 lb. weight capacity and our new comfort spring suspension for a smoother ride, yet it disassembles quickly and easily for transport. Standard features include a new ultra-bright LED angle-adjustable headlight, LED battery gauge, and more foot room than other scooters in its class. It's tight turning radius makes it ideal for indoor use, while it's ample ground clearance makes it great for outdoor use, as well.





## Wireless disassembly in 4 easy steps!

Item	GB117Z	GB147Z
Description	XLSHD	XLSHD
Medicare Code	Not Coded	Not Coded
Weight Capacity	325 lbs.	325 lbs.
Drive Wheels	Rear	Rear
Maximum Speed	4 mph	4 mph
Operating Range up to <sup>1</sup>	14 miles	14 miles
Under Deck Clearance	3"	3"
Turning Radius	37"	47"
Type Batteries	2-22 AH	2-22 AH
Colors	Red, Blue	Red, Blue
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brake	es Yes	Yes
Charger	Off-board	Off-board
Controller	50 Amp Dynamic R-Series	50 Amp Dynamic R-Series
Dimensions		
Length	42.5"	41.5"
Width	21.5"	21.5"
Height (ground to top of back of seat)	33"-35"	33"-35"
Ground to Top of Seat	21"-23"	21"-23"
Ground to Top of Deck	4.75"	4.75"
Deck to Top of Seat	16"-18"	16"-18"
Number of Seat Height P		3
Size of Increments	1"	1"
Front Axle to Rear Axle	33.5"	32"
Weight of Unit: (assembled) <sup>3</sup>	114 lbs.	122 lbs.
Front Half	35 lbs.	43 lbs.
Rear Half	28 lbs.	28 lbs.
Seat w/Arms	22 lbs.	22 lbs.
Batteries	29 lbs. pack	29 lbs. pack
Tires:	Solid	Solid
Front	8"	8"
Rear	8"	8"
Rear Anti-Tip Wheels	2"	2"
Standard Seat	Stadium Style	Stadium Style
Seat Back Height - no headrest	16"	16"
Width x Depth	18" x 16"	18" x 16"
Color	Black Vinyl	Black Vinyl
Optional Seat:	Stadium Style	Stadium Style
Back Height - no headrest	16"	16"
Width x Depth	20"x 18"	20"x18"
Color	Black Vinyl	Black Vinyl
Warranty:		
Frame	Lifetime Warranty	Lifetime Warranty
Drive Train	2 Year Warranty	2 Year Warranty
Electronics <sup>2</sup>	2 Year Warranty	2 Year Warranty
Batteries	1 Year Warranty	1 Year Warranty

<sup>1.</sup> Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain, and battery type.



Step 1: Pop off the seat!



Step 2: Pull out the battery pack!



Step 3: Fold down tiller and lock in place!



**Step 4:** Pull forward on the drivetrain release lever to separate the front and rear frame sections!



Disassembles into 5 pieces!

4 wheel model shown



A battery charging port is conveniently located on the tiller under the control panel.

## **Features**

- 325 lb. weight capacity
- Tiller mounted basket
- Angle-adjustable LED headlight
- 22 amp battery
- Comfort Spring System

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



<sup>&</sup>lt;sup>2.</sup> Electronics warranty excludes batteries.