

ALL NEW GP208

# Alanté SPORT

- *Low-maintenance Color Through Shroud Design!*
- *Larger, Lightweight Composite Footplate!*
- *Generous High Back Captain's Seat*



# AlantéSPORT



## FEATURES

- Low center of gravity
- Most stable unit in its class
- Excellent maneuverability
- 28” turning radius
- 300 lb weight capacity
- Upgraded user friendly motor connections
- Auto shut-off charger system
- Electromagnetic braking system
- Padded flip up arms standard
- 9” driving wheels
- Larger 6” rear caster wheels

